






Third tubular plate 3.5 as bulge

angle-stable; for screws \varnothing 3.5 mm; plate thickness 1.1 mm; stainless steel



Item no. 4.8082.06

Item no.	Item description	EAN
4.8082.05	one third tubular plate with conical thread 3.5; 5 holes; 61 x 9 mm; stainless steel	 4 047075 119387
4.8082.06	one third tubular plate with conical thread 3.5; 6 holes; 73 x 9 mm; stainless steel	 4 047075 119394
4.8082.07	one third tubular plate with conical thread 3.5; 7 holes; 85 x 9 mm; stainless steel	 4 047075 119400
4.8082.08	one third tubular plate with conical thread 3.5; 8 holes; 97 x 9 mm; stainless steel	 4 047075 119417
4.8082.09	one third tubular plate with conical thread 3.5; 9 holes; 109 x 9 mm; stainless steel	 4 047075 119424