

Set Clavicle Plates 3.5

angle-stable; titanium; set no. 19.478.00




Item no. 5.3502.34, right



Item no. 5.3542.1418, right



Item no. 5.3522.08, right



Set no. 19.478.00

 Internal hexagon

Item no.	Item description	Quant.
Implants - plates		
Angle-stable lateral clavicle plate 3.5; titanium		
5.3502.34	3/4 holes; 65 x 10 mm; lateral; right	1
5.3502.36	3/6 holes; 82 x 10 mm; lateral; right	1
5.3512.34	3/4 holes; 65 x 10 mm; lateral; left	1
5.3512.36	3/6 holes; 82 x 10 mm; lateral; left	1

Angle-stable clavicle shaft plates 3.5; titanium		
5.3522.08	8 holes; 99 x 12 mm; right	1
5.3522.10	10 holes; 124 x 12 mm; right	1
5.3532.08	8 holes; 99 x 12 mm; left	1
5.3532.10	10 holes; 124 x 12 mm; left	1

AC-joint hook plates with conical thread 3.5; titanium		
5.3542.1418	offset 14; clasp 18; right	1
5.3542.1818	offset 18; clasp 18; right	1
5.3552.1418	offset 14; clasp 18; left	1
5.3552.1818	offset 18; clasp 18; left	1

Instruments		
10.500.27A	drill guide with handle for small fragment; handle length 55 mm	1
2.904.05	drill bit for quick coupling; Ø 2.5 mm x 105 mm; spiral length 50 mm; double spiral	1
2.9406.25	hex screwdriver with handle for screws Ø 2.7 mm - Ø 4.0 mm; AF 2.5 mm; conical; PPSU	1
2.953.60	gauge with clasp for screws with conical head thread Ø 3.5 mm; Ø 4.0 mm; measuring range 60 mm	1
2.954.01	screw forceps; self holding	1
2.977.01	drill guide for angle stable screwing; length 45 mm; small fragment	2
10.520.30 *	double drill guide Ø 2.5 mm/ Ø 3.5 mm, spring bushing	

* optional

Item no.	Item description	Quant.
Implants - screws		
angle-stable cortical screw Ø 3.5 mm; fully threaded; self tapping; titanium		
 Internal hexagon		
3.133.10	length 10 mm	5
3.133.12	length 12 mm	5
3.133.14	length 14 mm	5
3.133.16	length 16 mm	5
3.133.18	length 18 mm	5
3.133.20	length 20 mm	5
3.133.22	length 22 mm	5
3.133.24	length 24 mm	5


cortical screw Ø 3.5 mm; fully threaded; self tapping; titanium		
 Internal hexagon		
3.132.10	length 10 mm	3
3.132.12	length 12 mm	3
3.132.14	length 14 mm	3
3.132.16	length 16 mm	3
3.132.18	length 18 mm	3
3.132.20	length 20 mm	3

Wire		
6.031.16	Kirschner wire with trocar point and round end; Ø 1.6 mm x 150 mm; stainless steel	5

Container		
19.477.00	perforated autoclavable container with inset for instruments and implants for clavicle plate; small fragment	1

Set Clavicle Plates 3.5

angle-stable; titanium; set no. 19.478.00



Set no. 19.478.00



Optionally applicable:

WiKo® anterior medial clavicle plate 3.5; titanium; TiOB®

- | | |
|-----------|---|
| 5.3802.06 | WiKo® anterior medial clavicle plate 3.5; 6 hole; length 80 mm |
| 5.3802.07 | WiKo® anterior medial clavicle plate 3.5; 7 hole; length 91 mm |
| 5.3802.08 | WiKo® anterior medial clavicle plate 3.5; 8 hole; length 103 mm |

WiKo® anterior lateral clavicle plate 3.5; titanium; TiOB®

- | | |
|-----------|---|
| 5.3812.08 | WiKo® anterior lateral clavicle plate 3.5; 8 hole; length 105 mm |
| 5.3812.09 | WiKo® anterior lateral clavicle plate 3.5; 9 hole; length 117 mm |
| 5.3812.10 | WiKo® anterior lateral clavicle plate 3.5; 10 hole; length 128 mm |



Königsee Implantate GmbH
OT Aschau | Am Sand 4
07426 Allendorf | Germany

Fon +49 36738 498 - 560
Fax +49 36738 498 - 569

Mail export@koenigsee-implantate.de

For further information please visit:
www.koenigsee-implantate.com

HeBo | Edition 02b | 2020-08
All pictured products are CE compliant.