


# Set DHS/DCS

Instruments; Set no. 19.120.00



Set no. 19.120.00  Internal hexagon

Item no.	Item description	Quant.
<b>Instruments</b>		
2.904.08	drill bit for quick coupling; Ø3.2x145; spiral length 60; double spiral	1
2.904.12	drill bit for quick coupling; Ø4.5x145; spiral length 60; double spiral	1
10.500.02	double drill guide; Ø3.2/Ø4.5	1
10.500.21	drill guide for neutral and compression position; Ø3.2	1
2.9406.35	hex screwdriver with handle for screws Ø4.5-Ø7.0; AF3.5; length 215	1
2.941.35	hex screwdriver inset f.quick coupling f.screw Ø4.5 until Ø7.0; wrench width 3.5	1
2.953.45	gauge with clasp for screws Ø4.5 and Ø6.5; measuring range 115	1
10.100.01	gauge for DHS/DCS	1
10.152.01	universal targeting device for dynamic hip screws	1
10.300.01	T-handle with quick coupling for DHS; big	1

Item no.	Item description	Quant.
<b>Instruments</b>		
10.310.01	wrench for insertion and extraction of hip and condylar screws	1
10.310.02	guiding shank	1
10.315.05	tap for hip and condylar screws	1
10.320.01	three-level step drill for dynamic hip screws; complete	1
10.330.01	impactor for dynamic hip and condylar screws	1
10.340.01	guide sleeve; short	1
10.340.02	guide sleeve; long	1
10.350.01	connecting screw; length 145; short; st.st.	1
10.350.02	connecting screw; length 245; long; st.st.	1

#### Optionally available:

10.325.02	cutter for dynamic condylar screws
10.130.01	targeting device for dynamic hip screw 130°
10.135.01	targeting device for dynamic hip screw 135°
10.140.01	targeting device for dynamic hip screw 140°
10.145.01	targeting device for dynamic hip screw 145°
10.150.01	targeting device for dynamic hip screw 150°
10.321.01	targeting device for dynamic condylar screws

<b>Container</b>		
19.123.00	perforated autoclavable container with inset for DHS/DCS-instruments	1