


# Fracture compression screw 3.0/4.0 st sd Cann.

for DENS screwing; titanium



Fracture compressing screw  
Item no. 3.123.32

Item no.	Item description	EAN
----------	------------------	-----

Fracture compressing screw Ø 3.0 mm / Ø 4.0 mm for DENS screwing, inner Ø 1.3 mm; titanium

3.123.20	fracture compression screw;Ø3/Ø4x20;thread 6;sd;st;lØ1.3;T	 4047075088515
3.123.22	fracture compression screw;Ø3/Ø4x22;thread 6;sd;st;lØ1.3;T	 4047075088652
3.123.24	fracture compression screw;Ø3/Ø4x24;thread 6;sd;st;lØ1.3;T	 4047075088539
3.123.26	fracture compression screw;Ø3/Ø4x26;thread 6;sd;st;lØ1.3;T	 4047075088546
3.123.28	fracture compression screw;Ø3/Ø4x28;thread 6;sd;st;lØ1.3;T	 4047075088553
3.123.30	fracture compression screw;Ø3/Ø4x30;thread 6;sd;st;lØ1.3;T	 4047075088454
3.123.32	fracture compression screw;Ø3/Ø4x32;thread 6;sd;st;lØ1.3;T	 4047075088805
3.123.34	fracture compression screw;Ø3/Ø4x34;thread 6;sd;st;lØ1.3;T	 4047075088669
3.123.36	fracture compression screw;Ø3/Ø4x36;thread 6;sd;st;lØ1.3;T	 4047075088676
3.123.38	fracture compression screw;Ø3/Ø4x38;thread 6;sd;st;lØ1.3;T	 4047075088683
3.123.40	fracture compression screw;Ø3/Ø4x40;thread 6;sd;st;lØ1.3;T	 4047075088690
3.123.42	fracture compression screw;Ø3/Ø4x42;thread 6;sd;st;lØ1.3;T	 4047075088706
3.123.44	fracture compression screw;Ø3/Ø4x44;thread 6;sd;st;lØ1.3;T	 4047075088560

All implants are available on request as non-sterile packed or sterile packed.