

Set KI-CVA-plate

for screws Ø 4.5/5.2 mm; titanium; Set no. 19.651.00



Item no. 3.031.12




Item no. 3.033.12




Item no. 5.555.1624



Item no. 5.555.1652



Set no. 19.651.00  six lobe

Item no.	Item description	Quant.
Implants – plates		
KI-CVA-plates for screws Ø 4.5 mm/5.2 mm; titanium		
5.555.1624	4 holes; 16 x 24 mm	1
5.555.1626	4 holes; 16 x 26 mm	1
5.555.1628	4 holes; 16 x 28 mm	1
5.555.1630	6 holes; 16 x 30 mm	1
5.555.1632	6 holes; 16 x 32 mm	1
5.555.1634	6 holes; 16 x 34 mm	1
5.555.1636	6 holes; 16 x 36 mm	1
5.555.1638	6 holes; 16 x 38 mm	1
5.555.1640	6 holes; 16 x 40 mm	1
5.555.1642	6 holes; 16 x 42 mm	1
5.555.1644	6 holes; 16 x 44 mm	1
5.555.1646	8 holes; 16 x 46 mm	1
5.555.1648	8 holes; 16 x 48 mm	1
5.555.1650	8 holes; 16 x 50 mm	1
5.555.1652	8 holes; 16 x 52 mm	1
5.555.1654	8 holes; 16 x 54 mm	1
5.555.1657	8 holes; 16 x 57 mm	1
5.555.1660	8 holes; 16 x 60 mm	1

Item no.	Item description	Quant.
Implants – screws		
Spine screws Ø 4.5 mm; monocortically; titanium		
3.031.12	length 12 mm	8
3.031.14	length 14 mm	8
3.031.16	length 16 mm	8
3.031.18	length 18 mm	8
3.031.20	length 20 mm	8
Spine screws for revision; Ø 5.2 mm; monocortically; titanium		
3.033.12	length 12 mm	4
3.033.14	length 14 mm	4
3.033.16	length 16 mm	4
3.033.18	length 18 mm	4
3.033.20	length 20 mm	4

Instruments		
10.268.51	cone point ball; ball Ø 6.0 mm	1
10.268.52	awl; Ø 2.5 mm; conical	1
10.268.53	screwdriver; six-lobe 10	1
10.268.54	roll-bending plier	1
10.268.06	screwdriver for distraction pins	1
10.268.07	fixation pin Ø 3.5 mm; wrench size 4.0 mm; length 10 mm	5


Container		
19.650.00	perforated autoclavicle container with insert for instruments/implants for KI-CVA-plates	1


Information:

In addition to the standard set KI-CVA-plates up to a length of 93 mm can be ordered for level 2.

Set KI-CVA-plate

for screws \varnothing 4.5/5.2 mm; titanium; Set no. 19.651.00



Set no. 19.651.00  six lobe

