

# Set L- and T-Tibia head plates lateral as

angle-stable; for screws Ø 6.0/4.5 mm; large fragment; titanium; Set no. 19.407.00



Item no. 5.773.46




Item no. 5.777.46



Item no. 5.774.46



Item no. 5.778.46



Set no. 19.407.00 Internal hexagon

Item no.	Item description	Quant.
<b>Implants – plates</b>		
<b>L-tibia head plate with conical thread 6.0/4.5; titanium</b>		
5.773.44	4 holes; 95 x 14 mm; left	1
5.773.46	6 holes; 127 x 14 mm; left	1
5.773.48	8 holes; 159 x 14 mm; left	1
5.773.410	10 holes; 191 x 14 mm; left	1
5.774.44	4 holes; 95 x 14 mm; right	1
5.774.46	6 holes; 127 x 14 mm; right	1
5.774.48	8 holes; 159 x 14 mm; right	1
5.774.410	10 holes; 191 x 14 mm; right	1

<b>T-tibia head plate with conical thread 6.0/4.5; titanium</b>		
5.777.44	4 holes; 95 x 14 mm; left	1
5.777.46	6 holes; 127 x 14 mm; left	1
5.777.48	8 holes; 159 x 14 mm; left	1
5.777.410	10 holes; 191 x 14 mm; left	1
5.778.44	4 holes; 95 x 14 mm; right	1
5.778.46	6 holes; 127 x 14 mm; right	1
5.778.48	8 holes; 159 x 14 mm; right	1
5.778.410	10 holes; 191 x 14 mm; right	1


<b>Implants – screws</b>		
<b>cancellous screw with conical head thread; Ø 6.0 mm; fully threaded; self tapping; titanium</b> Internal hexagon		
3.163.55	length 55 mm	4
3.163.60	length 60 mm	4
3.163.65	length 65 mm	4

Item no.	Item description	Quant.
<b>Implants – screws</b>		
<b>cancellous screw with conical head thread; Ø 6.0 mm; fully threaded; self tapping; titanium</b> Internal hexagon		
3.163.70	length 70 mm	4
3.163.75	length 75 mm	4
3.163.80	length 80 mm	4
3.163.85	length 85 mm	4
3.163.90	length 90 mm	4

<b>cortical screw with conical head thread; Ø 4.5 mm; fully threaded; self tapping; titanium</b> Internal hexagon		
3.545.18	length 18 mm	5
3.545.20	length 20 mm	5
3.545.22	length 22 mm	5
3.545.24	length 24 mm	5
3.545.26	length 26 mm	5
3.545.28	length 28 mm	5
3.545.30	length 30 mm	5
3.545.32	length 32 mm	5
3.545.34	length 34 mm	5
3.545.36	length 36 mm	5
3.545.38	length 38 mm	5
3.545.40	length 40 mm	5
3.545.42	length 42 mm	5
3.545.45	length 45 mm	5
3.545.50	length 50 mm	5
3.545.55	length 55 mm	5

# Set L- and T-Tibia head plates lateral as

angle-stable; for screws Ø 6.0/4.5 mm; large fragment; titanium; Set no. 19.407.00





Set no. 19.407.00  Internal hexagon

Item no.	Item description	Quant.
<b>Instruments</b>		
2.9406.35	hex screwdriver with handle for screws Ø 4.5 mm - Ø 7.0 mm, AF 3.5, length 215 mm	1
2.904.09	drill bit for quick coupling, Ø 3.2 x 195 mm, spiral length 60 mm, double spiral	1
2.977.02	drill guide for angle stable screwing, length 70 mm, big fragment	2
2.953.90	gauge with clasp for screws Ø 4.5 mm, Ø 6.0 mm, measuring range 90 mm	1
2.954.01	screw forceps, self holding	1

<b>Wires</b>		
6.031.18	Kirschner wire with trocar point and round end, Ø 1.8 x 150 mm; stainless steel	5

<b>Container</b>		
19.400.00/8	perforated autoclavable container with inset for instruments and proximal tibia head plates; large fragment	1

Item no.	Item description	Quant.
<b>Implants – screws</b>		
cortical screw Ø 4.5 mm; fully threaded; self tapping; titanium 		
3.152.18	length 18 mm	5
3.152.20	length 20 mm	5
3.152.22	length 22 mm	5
3.152.24	length 24 mm	5
3.152.26	length 26 mm	5
3.152.28	length 28 mm	5
3.152.30	length 30 mm	5
3.152.32	length 32 mm	5
3.152.34	length 34 mm	5
3.152.36	length 36 mm	5
3.152.38	length 38 mm	5
3.152.40	length 40 mm	5
3.152.42	length 42 mm	5
3.152.45	length 45 mm	5
3.152.50	length 50 mm	5