

Set Dist. radius plate 2.7 vas IH T / IMPLANTS sterile

variable angle-stable; for screws Ø 2.7 mm; with internal hexagon 2.5; titanium

not sterile packed Set no. 19.1163.00 / sterile packed Set no. 19.1163.00PS+



Item no. 5.9823.53
Item no. S5.9823.53




Item no. 5.9833.53
Item no. S5.9833.53



Item no. 5.9843.73
Item no. S5.9843.73



Item no. 5.9853.73
Item no. S5.9853.73



All implants mentioned are available individually, not sterile and sterile packed.

Implants in the set **not sterile packed** Set no. 19.1163.00

Implants in the set **sterile packed** Set no. 19.1163.00PS+

Internal hexagon
Set no. 19.1163.00
not sterile packed

Item no.	Item no. sterile	Item description	Quant.
Implants – plates			
Distal radius plate 2.7 vas narrow			
5.9823.53	S5.9823.53	Distal radius plate with conic thread Ø 2.7 mm; 5/3 holes; 51 mm x 10 mm; narrow; right; variable angle-stable; titanium	1
5.9823.54	S5.9823.54	Distal radius plate with conic thread Ø 2.7 mm; 5/4 holes; 60 mm x 10 mm; narrow; right; variable angle-stable; titanium	1
5.9833.53	S5.9833.53	Distal radius plate with conic thread Ø 2.7 mm; 5/3 holes; 51 mm x 10 mm; narrow; left; variable angle-stable; titanium	1
5.9833.54	S5.9833.54	Distal radius plate with conic thread Ø 2.7 mm; 5/4 holes; 60 mm x 10 mm; narrow; left; variable angle-stable; titanium	1
Distal radius plate 2.7 vas broad			
5.9843.73	S5.9843.73	Distal radius plate with conic thread Ø 2.7 mm; 7/3 holes; 51 mm x 10 mm; broad; right; variable angle-stable; titanium	1
5.9843.74	S5.9843.74	Distal radius plate with conic thread Ø 2.7 mm; 7/4 holes; 60 mm x 10 mm; broad; right; variable angle-stable; titanium	1
5.9853.73	S5.9853.73	Distal radius plate with conic thread Ø 2.7 mm; 7/3 holes; 51 mm x 10 mm; broad; left; variable angle-stable; titanium	1
5.9853.74	S5.9853.74	Distal radius plate with conic thread Ø 2.7 mm; 7/4 holes; 60 mm x 10 mm; broad; left; variable angle-stable; titanium	1

Implants – screws			
Cortical screw 2.7 as st 			
3.125.12	S3.125.12	Cortical screw; Ø 2.7 mm x 12 mm; angle-stable; fully threaded; self tapping; head Ø 5.0 mm; internal hexagon 2.5; titanium	5
3.125.14	S3.125.14	Cortical screw; Ø 2.7 mm x 14 mm; angle-stable; fully threaded; self tapping; head Ø 5.0 mm; internal hexagon 2.5; titanium	5
3.125.16	S3.125.16	Cortical screw; Ø 2.7 mm x 16 mm; angle-stable; fully threaded; self tapping; head Ø 5.0 mm; internal hexagon 2.5; titanium	5
3.125.18	S3.125.18	Cortical screw; Ø 2.7 mm x 18 mm; angle-stable; fully threaded; self tapping; head Ø 5.0 mm; internal hexagon 2.5; titanium	5
3.125.20	S3.125.20	Cortical screw; Ø 2.7 mm x 20 mm; angle-stable; fully threaded; self tapping; head Ø 5.0 mm; internal hexagon 2.5; titanium	5

Set Dist. radius plate 2.7 vas IH T / IMPLANTS sterile

variable angle-stable; for screws Ø 2.7 mm; with internal hexagon 2.5; titanium
not sterile packed Set no. 19.1163.00 / sterile packed Set no. 19.1163.00PS+

Item no.	Item no. sterile	Item description	Quant.
3.125.22	S3.125.22	Cortical screw; Ø 2.7 mm x 22 mm; angle-stable; fully threaded; self tapping; head Ø 5.0 mm; internal hexagon 2.5; titanium	5
3.125.24	S3.125.24	Cortical screw; Ø 2.7 mm x 24 mm; angle-stable; fully threaded; self tapping; head Ø 5.0 mm; internal hexagon 2.5; titanium	5
3.125.26	S3.125.26	Cortical screw; Ø 2.7 mm x 26 mm; angle-stable; fully threaded; self tapping; head Ø 5.0 mm; internal hexagon 2.5; titanium	5
3.125.28	S3.125.28	Cortical screw; Ø 2.7 mm x 28 mm; angle-stable; fully threaded; self tapping; head Ø 5.0 mm; internal hexagon 2.5; titanium	5

Cortical screw 2.7 st



3.120.12	S3.120.12	Cortical screw; Ø 2.7 mm x 12 mm; fully threaded; self tapping; internal hexagon 2.5; titanium	5
3.120.14	S3.120.14	Cortical screw; Ø 2.7 mm x 14 mm; fully threaded; self tapping; internal hexagon 2.5; titanium	5
3.120.16	S3.120.16	Cortical screw; Ø 2.7 mm x 16 mm; fully threaded; self tapping; internal hexagon 2.5; titanium	5
3.120.18	S3.120.18	Cortical screw; Ø 2.7 mm x 18 mm; fully threaded; self tapping; internal hexagon 2.5; titanium	5
3.120.20	S3.120.20	Cortical screw; Ø 2.7 mm x 20 mm; fully threaded; self tapping; internal hexagon 2.5; titanium	5

Instruments

2.904.03	S2.904.03	Drill bit for quick coupling Ø 2.0 mm; spiral 30 mm; 2-flute; length 105 mm; orange	1
2.9406.25	—	Hex screwdriver with handle for screws Ø 2.7 mm - 4.0 mm; wrench width 2.5 mm; conical	1
2.953.65	—	Gauge with clasp for screws Ø 2.7 mm; measuring range 60 mm	1
2.954.01	—	Screw forceps; self holding	1
2.977.19	—	Drill guide with handle for variable angle-stable screwing 2.7; length 42 mm	1
2.977.11/3	—	Drill guide for angle-stable screwing 2.7; length 44 mm; single thread with knurl	2

Wires

6.031.16	S6.031.16	Kirschner wire with trocar point and round end; Ø 1.6 mm x 150 mm; stainless steel	5
----------	-----------	--	---

Container

19.115.00	—	Perforated autoclavable container with inset for instruments/implants radius set 2.7	1
-----------	---	--	---


Optional

10.180.05	—	Aiming block for distal radius plate 5.9823.xx; right	—
10.180.07	—	Aiming block for distal radius plate 5.9843.xx; right	—
10.181.05	—	Aiming block for distal radius plate 5.9833.xx; left	—
10.181.07	—	Aiming block for distal radius plate 5.9853.xx; left	—

Set Dist. radius plate 2.7 vas IH T / IMPLANTS sterile

variable angle-stable; for screws Ø 2.7 mm; with internal hexagon 2.5; titanium
not sterile packed Set no. 19.1163.00 / sterile packed Set no. 19.1163.00PS+

Item no.	Item no. sterile	Item description	Quant.
Optional			
2.9345.25	—	Torque screwdriver handle; 2.2 newtonmeter	—
2.9355.35/20	—	Silicone handle with ratchet and AO coupling; length 165 mm; Ø 20 mm	—
2.941.25A	—	Hex screwdriver inset for quick coupling for screws Ø 2.7 mm – 4.0 mm; wrench width 2.5 mm; conical; length 100 mm	—



 Internal hexagon
Set no. 19.1163.00
not sterile packed

Set no. 19.1163.00PS+ applicable with: 19.027.00I MF 2.7 tools IH INSTRUMENTS
19.1163.00 Dist. radius plate 2.7 vas IH T SET