


# RoSA® Implant removal

Set no. 19.610.00




Set no. 19.610.00  Internal hexagon

Item no.	Item description	Quant.
Instruments		
2.9406.35/260	hex screwdriver with handle for screws Ø 4.5 mm - Ø 7.0 mm, AF 3.5 mm, length 260 mm	1
10.269.80	holding sleeve Ø 12 x 130 mm for RoSA®-connecting screw	1
10.269.63	Extraction rod for RoSA®-blade	1
10.269.631	Mass for extraction rod for RoSA®-blade	1
10.269.632	Screw hex 3.5 mm for mass for extraction rod RoSA®-blade	1
10.269.41	Screwdriver with T-handle for removal RoSA®-screw	1
10.269.50	threaded rod M4 x 305.5 mm for RoSA®, stainless steel	1
10.269.51	screwing aid for rod M4 for RoSA®	1
10.269.72	compression nut, M4 for RoSA®	1

Item no.	Item description	Quant.
Container		
19.609.00	perforated autoclavable container with inset for implant removal RoSA®	1

# RoSA® Implant removal

Set no. 19.610.00



Set no. 19.610.00  Internal hexagon