

Instruments RoSA®

Set no. 19.608.00



Set no. 19.608.00  Internal hexagon


Item no.	Item description	Quant.
Instruments		
Level 1		
10.305.06	T-handle with 3 points coupling; length 85 mm	1
10.269.122	angle guide 122° for RoSA®	1
10.269.129	angle guide 129° for RoSA®	1
10.269.136	angle guide 136° for RoSA®	1
10.269.143	angle guide 143° for RoSA®	1
10.269.10	wire drill guide for RoSA®	1
10.269.11	impactor for K-wires Ø 3.0 mm for RoSA®	1
10.421.05	extraction forceps with notch for wires	1
10.269.20	direct measuring gauge for RoSA®; length 260 mm; measuring range up to 170 mm	1
10.269.30	RoSA® quadruple reamer	1
10.269.40	insertion instrument for RoSA®-screw	1

Item no.	Item description	Quant.
Wires		
6.031.30/270	Kirschner wire with trocar point and round end; Ø 3.0 x 270 mm for RoSA®; stainless steel	4
6.031.30/400	Kirschner wire with trocar point and round end; Ø 3.0 x 400 mm for RoSA®; stainless steel	2
6.031.30/320	Kirschner wire with trocar point and round end; Ø 3.0 x 320 mm for RoSA®; stainless steel	2

Container		
19.607.00	perforated autoclavable container with inset for instruments RoSA®	1

Instruments RoSA®

Set no. 19.608.00




Set no. 19.608.00  Internal hexagon

Item no.	Item description	Quant.
Instruments		
Level 2		
10.269.50	threaded rod M4 x 305.5 mm for RoSA®; stainless steel	1
10.269.51	screwing aid for rod M4 for RoSA®	1
10.269.60	guide sleeve Ø 6.0 x 260 mm for RoSA®-blade	1
10.269.61	Fixing sleeve M4 for guide sleeve Ø 6.0 x 260 mm for RoSA®	2
10.269.62	impactor for RoSA®-blade	1
10.264.24	plastic hammer	1
10.269.70	compression instrument with guiding pin for RoSA®	1
10.269.7115	distance sleeve Ø10 x 15 mm for RoSA®	2
10.269.7110	distance sleeve Ø 10 x 10 mm for RoSA®	2
10.269.72	compression nut; M4 for RoSA®	1
10.269.80	holding sleeve Ø 12 x 130 mm for RoSA®-connecting screw	1
10.269.81	guiding sleeve for Redon tube CH8 for RoSA®	1
10.269.82	drill bit for quick coupling; Ø 3.2 x 305 mm; spiral length 55 mm; with marking	1
10.269.83	universal drill guide with scaling; Ø 3.2 mm for RoSA®	1


Item no.	Item description	Quant.
Instruments		
10.269.84	drill guide with scaling for angle stable screwing for drill Ø 3.2 mm	1
2.953.90	gauge with clasp for screws Ø 4.5 mm and Ø 6.0 mm; measuring range 90 mm	1
2.927.45	tap for quick coupling for cortical screws; Ø 4.5 x 195 mm; thread 70 mm	1
2.954.01	screw forceps; self holding	1
2.9406.35/260	hex screwdriver with handle for screws Ø 4.5 mm - Ø 7.0 mm; AF 3.5 mm; length 260 mm	1
2.941.35/195	hex screwdriver inset for quick coupling for screw Ø 4.5 mm - 7.0 mm; length 195 mm; AF 3.5 mm	1
optional:		
2.9771.45	drill guide for angle stable compression hole; big fragment	1
2.977.02	drill guide for angle stable screwing; length 70 mm; big fragment	1

Instruments RoSA®

Set no. 19.608.00




Set no. 19.608.00, upper level (1)  Internal hexagon



Set no. 19.608.00, lower level (2)  Internal hexagon

Instruments RoSA®

Set no. 19.608.00



Set no. 19.608.00  Internal hexagon