


Upgrade Set for TIS™-Screws

with Instruments; TiOB; Titanium; Set no. 19.681.00/1



Set no. 19.681.00/1

Item no.	Item description	Quant.
TIS™-screws; fully threaded		
Ø 7.5 mm, inside Ø 2.9 mm; self drilling; self tapping; cannulated; titanium		
3.530.135	length 135 mm	2
3.530.140	length 140 mm	2
3.530.145	length 145 mm	2
3.530.150	length 150 mm	2
3.530.155	length 155 mm	2
3.530.160	length 160 mm	2
3.530.165	length 165 mm	2

Item no.	Item description	Quant.
TIS™-screws; thread 32 mm		
Ø 7.5 mm, inside Ø 2.9 mm; self drilling; self tapping; back tapping; cannulated; titanium		
3.531.135	length 135 mm	2

Item no.	Item description	Quant.
Container		
19.680.00/1	perforated autoclavable container with inset for instruments and implants for upgrade-set TIS™-screws	1

Item no.	Item description	Quant.
Further instruments		
2.980.28/450	guide wire; Ø 2.8 x 450 mm; thread 5 mm; stainless steel	5

Information:

The name TIS™ derives from the German designation "Transiliosakralschraube", which translates into English as "transiliosacral screw". Products with the add-on "™" are trademarks of Königsee Implantate.

TiOB® is a surface treatment that brings optimized properties.