


Set VaWiKo® EPI plate 4.0 vas IH T INSTRUMENTS

for screws Ø 4.0 mm

non-sterile packaging Set no. 19.563.00I



 Internal hexagon
Set no. 19.563.00I
non-sterile packaging

Implants in sets, supplied **non-sterile** Set no. 19.563.00I

Item no.	Item no. sterile	Item description	Quant.
Instruments			
2.977.21	—	Drill guide with handle for variable angle-stable screwing 4.0 mm; for drills Ø 3.0 mm; length 68 mm	1
10.940.03	—	Drill guide for drill bush system 10.940.00; outside Ø 3.2 mm; internal Ø 1.9 mm	1
2.950.150	—	Gauge for cannulated screws Ø 2.7 mm – 4.0 mm; length 150 mm; measuring range 80 mm	1
2.906.40	S2.906.40	Pilot drill bit for quick coupling 3.0 mm; drill length 60 mm; cannulated; internal Ø 1.8 mm; scaled; length 100 mm	1
2.9436.40	—	Hex screwdriver with handle for Ø 2.7 mm – 4.0 mm; cannulated; internal Ø 1.7 mm; wrench width 2.5 mm; length 215 mm	1
2.977.16	—	Drill guide for angle-stable screwing 4.0 mm; for drills Ø 3.0 mm; length 68 mm	2
2.904.25	S.2.904.25	Drill bit for quick coupling 3.0 mm; spiral 60 mm; 2-flute; length 155 mm	1
2.977.14	—	Drill guide for variable angle-stable screwing 2.0 mm; for drills Ø 1.5 mm; length 45 mm	1
2.953.65	—	Gauge with clasp for screws Ø 2.7 mm; measuring range 60 mm	1


Wires			
2.982.16	—	Guide wire; Ø 1.6 mm x 150 mm; thread 10 mm; stainless steel	5
6.031.10	S.6.031.10	Kirschner wire with trocar point and round end; Ø 1.0 mm x 150 mm; stainless steel	5

Container			
19.562.00	—	Perforated autoclavable container with inset for VaWiKo® EPI plate and instruments	1

Set VaWiKo® EPI plate 4.0 vas IH T INSTRUMENTS

for screws Ø 4.0 mm

non-sterile packaging Set no. 19.563.00I



 Internal hexagon
Set no. 19.563.00I
non-sterile packaging

Information

The system name VaWiKo® is derived from the German designation „Variabel winkelstabile Kombination“, which translates into English as „variable angle-stable combination“. This system combines conventional plating technologies with angle-stable locking in variable angles.

VaWiKo® is a registered trademark of Königsee Implantate GmbH.



Königsee Implantate GmbH
OT Aschau | Am Sand 4
07426 Allendorf | Germany

Fon +49 36738 498-560
Fax +49 36738 498-569

export@koenigsee-implantate.de

For further information please visit:
www.koenigsee-implantate.com

KrMo | Edition 04 | 2026-02
All pictured products are CE compliant.