

# Set Instruments VaWiKo® EPI plate 4.0 vas IH

variable angle stable; for screws Ø 4.0 mm; Set no. 19.563.001



Set no. 19.563.001  Internal hexagon

Item no.	Item description	Quant.
<b>Instruments</b>		
2.977.17	drill guide with handle for variable angle stable screwing 4.0; length 68 mm	1
2.950.150	gauge for cannulated screws Ø 2.7 mm to Ø 4.0 mm; length 150 mm	1
2.906.40	Pilot Drill Bit qc 3.0; DL 60 mm; cannulated 1.8; scaled; length 100 mm	1
2.9436.40	hex screwdriver with handle for Ø 2.7 mm - 4.0 mm; cannulated; Inside Ø 1.7 mm; wrench size 2.5; length 215 mm	1
2.977.16	drill guide for angle stable screwing 4.0; length 68 mm	2
2.904.25	drill bit qc 3.0; Spiral 60 mm; 2-flute; length 155 mm	1

Item no.	Item description	Quant.
<b>Instruments</b>		
2.977.14	drill guide for variable angle stable screwing 2.0; L 45 mm	1
2.953.65	gauge with clasp for screws Ø 2.7 mm; measuring range 60 mm	1

<b>Wires</b>		
2.982.16	guide wire; Ø 1.6 x 150 mm; thread 10; stainless steel	5
6.031.10	Kirschner wire with trocar point and round end; Ø1.0 x 150mm; stainless steel	5

<b>Container</b>		
19.562.00	perforated autoclavable container with inset for VaWiKo® EPI-Plate and instruments	1