

# Set VaWiKo<sup>®</sup> proximal Humerus plate 4.0/3.5


variable angle-stable; for screws Ø 4.0/3.5 mm; titanium; set no. 19.3752.00



Item no. 5.07913.83






Item no. 5.07913.86



Set no. 19.3752.00 Internal hexagon


Item no.	Item description	Quant.
<b>Implants – plates</b>		
VaWiKo <sup>®</sup> proximal Humerus plate 4.0/3.5; vas; titanium		
5.07913.83	8/3 holes; length 77 mm	1
5.07913.84	8/4 holes; length 93 mm	1
5.07913.85	8/5 holes; length 109 mm	1
5.07913.86	8/6 holes; length 125 mm	1

<b>Implants – screws</b>		
Cortical screw Ø 3.5 mm; fully threaded; self tapping; titanium 		
3.132.22	length 22 mm	5
3.132.24	length 24 mm	5
3.132.26	length 26 mm	5
3.132.28	length 28 mm	5
3.132.30	length 30 mm	5
3.132.32	length 32 mm	5
3.132.34	length 34 mm	5
3.132.36	length 36 mm	5
3.132.38	length 38 mm	5
3.132.40	length 40 mm	5
3.132.42	length 42 mm	5

Item no.	Item description	Quant.
<b>Implants – screws</b>		
Cortical screw w. con. head thread; Ø 3.5 mm; fully thread.; self tapp.; titanium 		
3.133.22	length 22 mm	5
3.133.24	length 24 mm	5
3.133.26	length 26 mm	5
3.133.28	length 28 mm	5
3.133.30	length 30 mm	5
3.133.32	length 32 mm	5
3.133.34	length 34 mm	5
3.133.36	length 36 mm	5
3.133.38	length 38 mm	5
3.133.40	length 40 mm	5
3.133.42	length 42 mm	5
Cortical screw w. con. head thread; Ø 4.0 mm; self tapp.; titanium 		
3.148.34	length 34 mm	5
3.148.36	length 36 mm	5
3.148.38	length 38 mm	5
3.148.40	length 40 mm	5
3.148.42	length 42 mm	5
3.148.44	length 44 mm	5
3.148.46	length 46 mm	5
3.148.48	length 48 mm	5
3.148.50	length 50 mm	5
3.148.52	length 52 mm	5
3.148.54	length 54 mm	5
3.148.56	length 56 mm	5

# Set VaWiKo® proximal Humerus plate 4.0/3.5

variable angle-stable; for screws Ø 4.0/3.5 mm; titanium; set no. 19.3752.00



Set no. 19.3752.00  Internal hexagon

Item no.	Item description	Quant.
<b>Instruments</b>		
10.173.01	sleeve for guide block for proximal small fragment humerus plates	2
10.179.15	aiming block f. prox. humerus pl. 5.07913.8X	1
2.904.22	drill bit for quick coupling, Ø 2.5 x 155 mm, spiral length 60 mm, double spiral	1
2.904.25	drill bit for quick coupling, Ø 3.0 x 155 mm, spiral length 60 mm, double spiral	1
2.9405.25	hex screwdriver for Ø 2.7 mm - Ø 4.0 mm; wrench size 2.5 mm; conical; length 215 mm	1
2.953.70	gauge with clasp for screws Ø 3.5 mm, Ø 4.0 mm, measuring range 70 mm	1
2.954.01	screw forceps, self holding	1
2.977.04	drill guide for angle stable screwing, length 68 mm, small fragment	1
2.977.09	drill guide w.handle for variable angle stable screwing, small fragment, length 68 mm	1
2.977.16	drill guide for angle stable screwing 4.0, length 68 mm	1
2.977.17	drill guide with handle for variable angle stable screwing 4.0, length 68 mm	1

Item no.	Item description	Quant.
<b>Wire</b>		
6.031.16	Kirschner wire with trocar point and round end, Ø 1.6 mm x 150 mm, stainless steel	5

<b>Container</b>		
19.3742.00	perforated autoclavable container with inset for VaWiKo® proximal humerus plates 5.07913.8X; Ø 4.0 mm; Ø 3.5 mm	1

## Information

The system name VaWiKo® is derived from the German designation „Variabel winkelstabile Kombination“, which translates into English as „variable angle-stable combination“. This system combines conventional plating technologies with angle-stable locking in variable angles.

VaWiKo® is a registered trademark of Königsee Implantate GmbH.



Königsee Implantate GmbH  
OT Aschau | Am Sand 4  
07426 Allendorf | Germany

Fon +49 36738 498-560  
Fax +49 36738 498-569

Mail [export@koenigsee-implantate.de](mailto:export@koenigsee-implantate.de)

For further information please visit:  
[www.koenigsee-implantate.com](http://www.koenigsee-implantate.com)

HeBo | Edition 02 | 2019-10  
All pictured products are CE compliant.