


# Set Cable osteosynthesis Ø1.0


by Prof. Labitzke; stainless steel; Set no. 19.518.00



Item no. 6.409.10/150

Item no. 6.409.10/500

Item no. 6.409.P10



Set no. 19.518.00

Item no.	Item description	Quant.
<b>Implants – cable</b>		
6.409.10/150	cable with ear for cable osteosynthesis, Ø 1.0 x 150 mm, stainless steel	10
6.409.10/500	cable for cable osteosynthesis, Ø 1.0 x 500 mm, laser welded on both sides, stainless steel	10

<b>Implants – ferrule</b>		
6.409.P10	ferrule for cable Ø 1.0 mm, stainless steel	10

<b>Instruments</b>		
10.265.01	awl, bent	1
10.265.02	awl, straight	1
10.265.05	cable cutter	1
10.265.06	rope retainer	1
10.265.10	ferrule forceps for cable Ø 1.0 mm	1

Item no.	Item description	Quant.
<b>Container</b>		
19.517.00	perforated autoclavable container with inset for instruments for cable plate osteosynthesis	1

