


# Set CPL cable osteosynthesis Ø1.5


Cable-Plate-System Labitzke; stainless steel; Set no. 19.518.00/1



Item no. 6.409.15/500



Item no. 6.409.P15



Set no. 19.518.00/1

| Item no.                | Item description  | Quant. |
|-------------------------|---|--------|
| <b>Implants – cable</b> |   |        |
| 6.409.15/500            | cable for cable osteosynthesis; Ø 1.5 x 500 mm; laser welded on both sides; stainless steel | 10     |

| <b>Implants – ferrule</b> |   |    |
|---------------------------|---|----|
| 6.409.P15                 | ferrule for cable Ø 1.5 mm; stainless steel | 10 |

| <b>Instruments</b> |   |   |
|--------------------|---|---|
| 10.2656.03/50      | loop instrument for cable; 120° angled; Ø 50 mm | 1 |
| 10.2656.04/50      | loop instrument for cable; straight; Ø 50 mm    | 1 |
| 10.265.05          | cable cutter                                    | 1 |
| 10.265.06          | rope retainer                                   | 1 |
| 10.265.15          | ferrule forceps for cable Ø 1.5 mm              | 1 |


| Item no.         | Item description  | Quant. |
|------------------|---|--------|
| <b>Container</b> |   |        |
| 19.517.00        | perforated autoclavable container with inset for instruments for cable plate osteosynthesis 1.5 | 1      |

---

# Set CPL cable osteosynthesis Ø1.5

Cable-Plate-System Labitzke; stainless steel; Set no. 19.518.00/1

---



Set no. 19.518.00/1