

Set Dynamic Epiphyseal Telescopic Screw (DET)


Stainless Steel; Set no. 19.501.00




Item no. 6.285.80




Item no. 6.286.80



Item no. 19.501.00

Item no.	Item description	Quant.
Implants - screws		
DET-screw Ø 7.0 mm thread 7.0 mm; inside Ø 2.2 mm selftapping; backtapping; stainless steel		
6.285.40	length 40 mm	3
6.285.50	length 50 mm	3
6.285.60	length 60 mm	3
6.285.70	length 70 mm	3
6.285.80	length 80 mm	3

DET-screw sleeve Ø 7.5 mm inside Ø 5.6 mm; stainless steel		
6.286.40	length 40 mm	3
6.286.50	length 50 mm	3
6.286.60	length 60 mm	3
6.286.70	length 70 mm	3
6.286.80	length 80 mm	3


Wire		
2.988.20/310	guide wire; Ø 2.0 mm x 310 mm; thread 15 mm; stainless steel	5

Container		
19.500.00	perforated autoclavable container with inset for instruments and DET-screws	1

Item no.	Item description	Quant.
Instruments		
2.9436.35/275	hex screwdriver with handle for Ø 4.5 mm - Ø 7.0 mm; cannulated; inside Ø 2.2 mm; wrench size 3.5 mm; length 275 mm	1
2.9476.02	screwdriver for extraction for DET-screw; length 275 mm	1
2.908.76	drill bit with triangular shank; Ø 7.6 mm x 285 mm; spiral length 100 mm; cannulated; inside Ø 2.2 mm; triple spiral	1
10.261.19	drill bit with triangular shank; Ø 4.8 mm x 320 mm; triple spiral; spiral length 55 mm; cannulated; inside Ø 2.2 mm	1
10.261.01	drill guide; Ø 8.0 x 150 mm; green	1
10.261.03	drill guide; Ø 4.8 x 160 mm; blue	1
10.261.02	drill guide; Ø 2.1 x 160 mm; silver	1
2.950.02	gauge for DET-screw; length 145 mm; measuring range 120 mm	1
2.928.01	instrument for holding of wires for DET-screw; Ø 2.0 mm; length 347 mm	1
2.954.01	screw forceps; self holding	1

Set Dynamic Epiphyseal Telescopic Screw (DET)

Stainless Steel; Set no. 19.501.00



Item no. 19.501.00