

# Cancellous screw

## angle-stable; titanium



Item no. 3.147.40



Item no. 3.163.45

Item no.	Item description	EAN
<b>Cancellous screw Ø 4.0 mm; angle-stable; titanium</b>		
3.147.10	Cancellous screw with conical head thread; Ø 4.0 x 10 mm; fully threaded; self tapping; titanium	4 047075 091225
3.147.12	Cancellous screw with conical head thread; Ø 4.0 x 12 mm; fully threaded; self tapping; titanium	4 047075 091355
3.147.14	Cancellous screw with conical head thread; Ø 4.0 x 14 mm; fully threaded; self tapping; titanium	4 047075 091249
3.147.16	Cancellous screw with conical head thread; Ø 4.0 x 16 mm; fully threaded; self tapping; titanium	4 047075 091256
3.147.18	Cancellous screw with conical head thread; Ø 4.0 x 18 mm; fully threaded; self tapping; titanium	4 047075 091263
3.147.20	Cancellous screw with conical head thread; Ø 4.0 x 20 mm; fully threaded; self tapping; titanium	4 047075 091140
3.147.22	Cancellous screw with conical head thread; Ø 4.0 x 22 mm; fully threaded; self tapping; titanium	4 047075 091508
3.147.24	Cancellous screw with conical head thread; Ø 4.0 x 24 mm; fully threaded; self tapping; titanium	4 047075 091362
3.147.26	Cancellous screw with conical head thread; Ø 4.0 x 26 mm; fully threaded; self tapping; titanium	4 047075 091379
3.147.28	Cancellous screw with conical head thread; Ø 4.0 x 28 mm; fully threaded; self tapping; titanium	4 047075 091386
3.147.30	Cancellous screw with conical head thread; Ø 4.0 x 30 mm; fully threaded; self tapping; titanium	4 047075 091270
3.147.32	Cancellous screw with conical head thread; Ø 4.0 x 32 mm; fully threaded; self tapping; titanium	4 047075 091287
3.147.34	Cancellous screw with conical head thread; Ø 4.0 x 34 mm; fully threaded; self tapping; titanium	4 047075 091416
3.147.35	Cancellous screw with conical head thread; Ø 4.0 x 35 mm; fully threaded; self tapping; titanium	4 047075 091423
3.147.36	Cancellous screw with conical head thread; Ø 4.0 x 36 mm; fully threaded; self tapping; titanium	4 047075 091294
3.147.38	Cancellous screw with conical head thread; Ø 4.0 x 38 mm; fully threaded; self tapping; titanium	4 047075 091300
3.147.40	Cancellous screw with conical head thread; Ø 4.0 x 40 mm; fully threaded; self tapping; titanium	4 047075 091317

# Cancellous screw

## angle-stable; titanium



Item no. 3.147.40



Item no. 3.163.45

Item no.	Item description	EAN
<b>Cancellous screw Ø 4.0 mm; angle-stable; titanium</b>		
3.147.42	Cancellous screw with conical head thread; Ø 4x42 mm; fully threaded; self tapping; titanium	4 047075 091539
3.147.44	Cancellous screw with conical head thread; Ø 4x44 mm; fully threaded; self tapping; titanium	4 047075 091546
3.147.45	Cancellous screw with conical head thread; Ø 4x45 mm; fully threaded; self tapping; titanium	4 047075 091553
3.147.46	Cancellous screw with conical head thread; Ø 4x46 mm; fully threaded; self tapping; titanium	4 047075 091393
3.147.48	Cancellous screw with conical head thread; Ø 4x48 mm; fully threaded; self tapping; titanium	4 047075 091409
3.147.50	Cancellous screw with conical head thread; Ø 4x50 mm; fully threaded; self tapping; titanium	4 047075 091577
3.147.52	Cancellous screw with conical head thread; Ø 4x52 mm; fully threaded; self tapping; titanium	4 047075 091430
3.147.54	Cancellous screw with conical head thread; Ø 4x54 mm; fully threaded; self tapping; titanium	4 047075 091447
3.147.55	Cancellous screw with conical head thread; Ø 4x55 mm; fully threaded; self tapping; titanium	4 047075 091454
3.147.56	Cancellous screw with conical head thread; Ø 4x56 mm; fully threaded; self tapping; titanium	4 047075 091461
3.147.58	Cancellous screw with conical head thread; Ø 4x58 mm; fully threaded; self tapping; titanium	4 047075 091478
3.147.60	Cancellous screw with conical head thread; Ø 4x60 mm; fully threaded; self tapping; titanium	4 047075 091485
3.147.65	Cancellous screw with conical head thread; Ø 4x65 mm; fully threaded; self tapping; titanium	4 047075 116041
3.147.70	Cancellous screw with conical head thread; Ø 4x70 mm; fully threaded; self tapping; titanium	4 047075 091348
3.147.75	Cancellous screw with conical head thread; Ø 4x75 mm; fully threaded; self tapping; titanium	4 047075 116058
3.147.80	Cancellous screw with conical head thread; Ø 4x80 mm; fully threaded; self tapping; titanium	4 047075 116065

# Cancellous screw


















## angle-stable; titanium



Item no. 3.147.40



Item no. 3.163.45

Item no.	Item description	EAN
<b>Cancellous screw; Ø 6.0 mm; angle-stable; titanium</b>		
3.163.24	Cancellous screw with conical head thread; Ø 6.0 x 24 mm; fully threaded; self tapping; titanium	 4 047075 093472
3.163.26	Cancellous screw with conical head thread; Ø 6.0 x 26 mm; fully threaded; self tapping; titanium	 4 047075 093359
3.163.28	Cancellous screw with conical head thread; Ø 6.0 x 28 mm; fully threaded; self tapping; titanium	 4 047075 093366
3.163.30	Cancellous screw with conical head thread; Ø 6.0 x 30 mm; fully threaded; self tapping; titanium	 4 047075 093496
3.163.32	Cancellous screw with conical head thread; Ø 6.0 x 32 mm; fully threaded; self tapping; titanium	 4 047075 093373
3.163.34	Cancellous screw with conical head thread; Ø 6.0 x 34 mm; fully threaded; self tapping; titanium	 4 047075 093588
3.163.35	Cancellous screw with conical head thread; Ø 6.0 x 35 mm; fully threaded; self tapping; titanium	 4 047075 093595
3.163.36	Cancellous screw with conical head thread; Ø 6.0 x 36 mm; fully threaded; self tapping; titanium	 4 047075 093601
3.163.38	Cancellous screw with conical head thread; Ø 6.0 x 38 mm; fully threaded; self tapping; titanium	 4 047075 093489
3.163.40	Cancellous screw with conical head thread; Ø 6.0 x 40 mm; fully threaded; self tapping; titanium	 4 047075 093625
3.163.42	Cancellous screw with conical head thread; Ø 6.0 x 42 mm; fully threaded; self tapping; titanium	 4 047075 093502
3.163.44	Cancellous screw with conical head thread; Ø 6.0 x 44 mm; fully threaded; self tapping; titanium	 4 047075 093519
3.163.45	Cancellous screw with conical head thread; Ø 6.0 x 45 mm; fully threaded; self tapping; titanium	 4 047075 093526
3.163.46	Cancellous screw with conical head thread; Ø 6.0 x 46 mm; fully threaded; self tapping; titanium	 4 047075 093533
3.163.48	Cancellous screw with conical head thread; Ø 6.0 x 48 mm; fully threaded; self tapping; titanium	 4 047075 093403
3.163.50	Cancellous screw with conical head thread; Ø 6.0 x 50 mm; fully threaded; self tapping; titanium	 4 047075 093731
3.163.52	Cancellous screw with conical head thread; Ø 6.0 x 52 mm; fully threaded; self tapping; titanium	 4 047075 093632

# Cancellous screw

## angle-stable; titanium



Item no. 3.147.40



Item no. 3.163.45

Item no.	Item description	EAN
<b>Cancellous screw; Ø 6.0 mm; angle-stable; titanium</b>		
3.163.54	Cancellous screw with conical head thread; Ø 6.0 x 54 mm; fully threaded; self tapping; titanium	4 047075 093649
3.163.55	Cancellous screw with conical head thread; Ø 6.0 x 55 mm; fully threaded; self tapping; titanium	4 047075 093656
3.163.56	Cancellous screw with conical head thread; Ø 6.0 x 56 mm; fully threaded; self tapping; titanium	4 047075 093663
3.163.58	Cancellous screw with conical head thread; Ø 6.0 x 58 mm; fully threaded; self tapping; titanium	4 047075 093540
3.163.60	Cancellous screw with conical head thread; Ø 6.0 x 60 mm; fully threaded; self tapping; titanium	4 047075 093687
3.163.62	Cancellous screw with conical head thread; Ø 6.0 x 62 mm; fully threaded; self tapping; titanium	4 047075 119615
3.163.64	Cancellous screw with conical head thread; Ø 6.0 x 64 mm; fully threaded; self tapping; titanium	4 047075 119622
3.163.65	Cancellous screw with conical head thread; Ø 6.0 x 65 mm; fully threaded; self tapping; titanium	4 047075 093557
3.163.66	Cancellous screw with conical head thread; Ø 6.0 x 66 mm; fully threaded; self tapping; titanium	4 047075 119639
3.163.68	Cancellous screw with conical head thread; Ø 6.0 x 68 mm; fully threaded; self tapping; titanium	4 047075 119646
3.163.70	Cancellous screw with conical head thread; Ø 6.0 x 70 mm; fully threaded; self tapping; titanium	4 047075 093564
3.163.72	Cancellous screw with conical head thread; Ø 6.0 x 72 mm; fully threaded; self tapping; titanium	4 047075 126361
3.163.74	Cancellous screw with conical head thread; Ø 6.0 x 74 mm; fully threaded; self tapping; titanium	4 047075 119653
3.163.75	Cancellous screw with conical head thread; Ø 6.0 x 75 mm; fully threaded; self tapping; titanium	4 047075 059515
3.163.76	Cancellous screw with conical head thread; Ø 6.0 x 76 mm; fully threaded; self tapping; titanium	4 047075 128808
3.163.78	Cancellous screw with conical head thread; Ø 6.0 x 78 mm; fully threaded; self tapping; titanium	4 047075 128716
3.163.80	Cancellous screw with conical head thread; Ø 6.0 x 80 mm; fully threaded; self tapping; titanium	4 047075 093571

# Cancellous screw

angle-stable; titanium



Item no. 3.147.40



Item no. 3.163.45

Item no.	Item description	EAN
<b>Cancellous screw; Ø 6.0 mm; angle-stable; titanium</b>		
3.163.82	Cancellous screw with conical head thread; Ø 6.0 x 82 mm; fully threaded; self tapping; titanium	4 <sup>0</sup> 047075 <sup>0</sup> 128815
3.163.84	Cancellous screw with conical head thread; Ø 6.0 x 84 mm; fully threaded; self tapping; titanium	4 <sup>0</sup> 047075 <sup>0</sup> 128822
3.163.85	Cancellous screw with conical head thread; Ø 6.0 x 85 mm; fully threaded; self tapping; titanium	4 <sup>0</sup> 047075 <sup>0</sup> 059522
3.163.86	Cancellous screw with conical head thread; Ø 6.0 x 86 mm; fully threaded; self tapping; titanium	4 <sup>0</sup> 047075 <sup>0</sup> 128839
3.163.90	Cancellous screw with conical head thread; Ø 6.0 x 90 mm; fully threaded; self tapping; titanium	4 <sup>0</sup> 047075 <sup>0</sup> 059300