



Set Cannulated screws 7.0 st IH INSTRUMENTS BASIC

Set no. 19.223.001



Set no. 19.223.001  Internal hexagon

Item no.	Item description	Quant.
Instruments		
10.970.00	percutaneous drill bush system for cannulated screws Ø 7.0 mm	1
2.908.45	Drill bit ts 4.5; Spiral length 40 mm; 3-flute; length 230 mm; Cann. 2.1; scaled	1
2.907.70	countersink; Ø 8.0 x 160 mm; cutter ln 19.5 mm; cannulated; inside Ø 2.1 mm; for cannulated screws Ø 6.5 mm	1
10.306.06	T-handle with quick coupling; length 93.5 mm; cannulated; inside Ø 2.5 mm	1
2.942.02	hex screwdriver with T-handle for screw Ø 4.5 mm until Ø 7.0 mm; cannulated; inside Ø 2.2 mm	1
2.954.65A	holding sleeve for screws Ø 6.5 mm and Ø 7.0 mm	1
2.952.230	gauge for cannulated screws Ø 6.5 mm and Ø 7.0 mm; length 230 mm	1
2.960.03	parallel targeting device for guide wire Ø 2.0 mm	1

Item no.	Item description	Quant.
Implants – washers		
3.200.04	washer outside Ø 16 mm / inside Ø 6.7 mm; titanium	5

Wires		
2.979.20	guide wire, Ø 2.0 x 230 mm; thread 5; stainless steel	5

Container		
19.224.00	perforated autoclavable container with inset for instruments/cannulated screws Ø 7.0 mm	1

combinable with:

19.223.00	Cannul. screws 7.0 st sd IH T SET
19.223.00PS+	Cannul. screws 7.0 st sd IH T SET SCREWS sterile

All implants are available on request as non-sterile packed or sterile packed.

