


Set Instruments Cannulated Screws 2.7 IH

for cannulated screws \varnothing 2.7 mm; Set no. 19.027.201



Set no. 19.027.201  Internal hexagon

Item no.	Item description	Quant.
Instruments		
10.927.00	percutaneous drill bush system for cannulated screws \varnothing 2.7	1
2.906.27	Pilot Drill Bit qc2.2;DL60;cann. 1.5;scaled;L100	1
2.907.27	countersink, \varnothing 5x160,cannulated,inside \varnothing 1.5,for cannulated screws \varnothing 2.7	1
10.306.06	T-handle with quick coupling,length 93.5,cannulated,inside \varnothing 2.5	1
2.9436.40	hex screwdriver with handle for \varnothing 2.7-4;cannulated;l \varnothing 1.7;wrench size 2.5;L 215mm	1
2.954.27A	holding sleeve for screws \varnothing 2.7 to \varnothing 3.5	1
2.950.150	gauge for cannulated screws \varnothing 2.7 to \varnothing 4.0,length 150	1

Item no.	Item description	Quant.
Wires		
2.978.14	guide wire, \varnothing 1.4x150,thread 5,st.st.	5


Item no.	Item description	Quant.
Container		
19.040.01	Perforated autoclavable container for instruments can. screws 2.7-4.5	1

Set Instruments Cannulated Screws 2.7 IH

for cannulated screws \varnothing 2.7 mm; Set no. 19.027.201

Can be combined with:

- 19.170.00 Cann. screws 2.7 st sd IH T SET
19.170.00PS+ Cann. screws 2.7 st sd IH T SCREWS sterile



Set no. 19.027.201  Internal hexagon