


# Set Instruments Cannulated Screws 4.0 IH

for cannulated screws  $\varnothing$  4.0 mm; Set no. 19.040.201



Set no. 19.040.201  Internal hexagon

Item no.	Item description	Quant.
Instruments		
10.940.00	percutaneous drill bush system for cannulated screws $\varnothing$ 4	1
2.906.40	Pilot Drill Bit qc3.0;DL60;cann. 1.8;scaled;L100	1
2.907.40	countersink, $\varnothing$ 6x160,cannulated,inside $\varnothing$ 1.8,for cannulated screws $\varnothing$ 4	1
10.306.06	T-handle with quick coupling,length 93.5,cannulated,inside $\varnothing$ 2.5	1
2.9436.40	hex screwdriver with handle for $\varnothing$ 2.7-4;cannulated;l $\varnothing$ 1.7;wrench size 2.5;l 215mm	1
2.954.35A	holding sleeve for screws $\varnothing$ 3.5 and $\varnothing$ 4.0	1
2.950.150	gauge for cannulated screws $\varnothing$ 2.7 to $\varnothing$ 4.0,length 150	1

Item no.	Item description	Quant.
Wires		
2.978.16	guide wire, $\varnothing$ 1.6x150,thread 5,st.st.	5


Container		
19.040.01	Perforated autoclavable container for instruments can. screws 2.7-4.5	1

# Set Instruments Cannulated Screws 4.0 IH

for cannulated screws  $\varnothing$  4.0 mm; Set no. 19.040.201

Can be combined with:

- 19.190.00 Cann. screws 4.0 st sd IH T SET  
19.190.00PS+ Cann. screws 4.0 st sd IH T SCREWS sterile



Set no. 19.040.201  Internal hexagon