

Set Instruments Cannulated Screws 4.5 IH

for cannulated screws \varnothing 4.5 mm; Set no. 19.045.201



Set no. 19.045.201  Internal hexagon

Item no.	Item description	Quant.
Instruments		
10.945.00	percutaneous drill bush system for cannulated screws \varnothing 4.5	1
2.905.32	Drill bit qc3.2; Spiral60; 4-flute; L160; Cann.1.75	1
2.905.45	Drill bit qc4.5; Spiral60; 4-flute; L160; Cann.1.8	1
2.907.45	countersink, \varnothing 6.5x160, cannulated, inside \varnothing 1.75, for cannulated screws \varnothing 4.5	1
10.306.06	T-handle with quick coupling, length 93.5, cannulated, inside \varnothing 2.5	1
2.9436.35	hex screwdriver with handle for \varnothing 4.5-7; cannulated; I \varnothing 2.2; wrench size 3.5; L 215mm	1
2.954.45A	holding sleeve for screws \varnothing 4.5	1
2.951.150	gauge for cannulated screws \varnothing 4.5, length 150	1

Item no.	Item description	Quant.
Wires		
2.978.16	guide wire, \varnothing 1.6x150, thread 5, st.st.	5


Item no.	Item description	Quant.
Container		
19.040.01	Perforated autoclavable container for instruments can. screws 2.7-4.5	1

Set Instruments Cannulated Screws 4.5 IH

for cannulated screws \varnothing 4.5 mm; Set no. 19.045.20I

Can be combined with:

19.203.00	Cann. screws 4.5 st sd IH T SET
19.203.00PS+	Cann. screws 4.5 st sd IH T SCREWS sterile
19.200.00	Cann. screws 4.5 st IH T SET
19.200.00PS+	Cann. screws 4.5 st IH T SCREWS sterile



Set no.. 19.045.20I  Internal hexagon