

Set Cannulated screws 3.5 st sd


self tapping; self drilling; titanium; Set no. 19.184.00



Item no. 3.504.30



Item no. 3.604.30



Set no. 19.184.00  Internal hexagon

Item no.	Item description	Quant.
Implants – screws		
Partial thread - Ø 3.5 mm; self drilling; back tapping; self tapping; titanium		
3.504.10	length 10 mm; thread 5	2
3.504.12	length 12 mm; thread 5	2
3.504.14	length 14 mm; thread 5	2
3.504.16	length 16 mm; thread 6	2
3.504.18	length 18 mm; thread 7	2
3.504.20	length 20 mm; thread 8	2
3.504.22	length 22 mm; thread 9	2
3.504.24	length 24 mm; thread 10	2
3.504.26	length 26 mm; thread 12	2
3.504.28	length 28 mm; thread 14	2
3.504.30	length 30 mm; thread 14	2
3.504.32	length 32 mm; thread 14	2
3.504.34	length 34 mm; thread 14	2
3.504.36	length 36 mm; thread 14	2
3.504.38	length 38 mm; thread 14	2
3.504.40	length 40 mm; thread 14	2
3.504.42	length 42 mm; thread 15	2
3.504.44	length 44 mm; thread 15	2
3.504.46	length 46 mm; thread 15	2
3.504.48	length 48 mm; thread 15	2
3.504.50	length 50 mm; thread 15	2
3.504.52	length 52 mm; thread 16	2
3.504.54	length 54 mm; thread 16	2
3.504.56	length 56 mm; thread 16	2
3.504.58	length 58 mm; thread 16	2
3.504.60	length 60 mm; thread 16	2
3.504.62	length 62 mm; thread 17	2
3.504.64	length 64 mm; thread 17	2
3.504.66	length 66 mm; thread 17	2
3.504.68	length 68 mm; thread 17	2
3.504.70	length 70 mm; thread 17	2

Item no.	Item description	Quant.
Implants – screws		
Full thread - Ø 3.5 mm; self drilling; back tapping; self tapping; titanium		
3.604.10	length 10 mm; thread 5	2
3.604.12	length 12 mm; thread 5	2
3.604.14	length 14 mm; thread 5	2
3.604.16	length 16 mm; thread 6	2
3.604.18	length 18 mm; thread 7	2
3.604.20	length 20 mm; thread 8	2
3.604.22	length 22 mm; thread 9	2
3.604.24	length 24 mm; thread 10	2
3.604.26	length 26 mm; thread 12	2
3.604.28	length 28 mm; thread 14	2
3.604.30	length 30 mm; thread 14	2
3.604.32	length 32 mm; thread 14	2
3.604.34	length 34 mm; thread 14	2
3.604.36	length 36 mm; thread 14	2
3.604.38	length 38 mm; thread 14	2
3.604.40	length 40 mm; thread 14	2
3.604.42	length 42 mm; thread 15	2
3.604.44	length 44 mm; thread 15	2
3.604.46	length 46 mm; thread 15	2
3.604.48	length 48 mm; thread 15	2
3.604.50	length 50 mm; thread 15	2
3.604.52	length 52 mm; thread 16	2
3.604.54	length 54 mm; thread 16	2
3.604.56	length 56 mm; thread 16	2
3.604.58	length 58 mm; thread 16	2
3.604.60	length 60 mm; thread 16	2
3.604.62	length 62 mm; thread 17	2
3.604.64	length 64 mm; thread 17	2
3.604.66	length 66 mm; thread 17	2
3.604.68	length 68 mm; thread 17	2
3.604.70	length 70 mm; thread 17	2

Set Cannulated screws 3.5 st sd


self tapping; self drilling; titanium; Set no. 19.184.00

Item no.	Item description	Quant.
Instruments		
2.950.150	gauge for cannulated screws Ø2.7 to Ø4.0,length 150	1
10.306.06	T-handle with quick coupling,length 93.5,cannulated,inside Ø2.5	1
10.935.00	percutaneous drill bush system for cannulated screws Ø3.5	1
20.502.01	holding clamp,small	1
2.905.27	drill bit f.quick coupling,Ø2.7x160,spiral ln 30,quadr.spiral,cann.,inside Ø1.35	1
2.905.35	drill bit f.quick coupling,Ø3.5x160,spiral ln 30,quadr.spiral,cann.,inside Ø1.35	1
2.907.35	countersink,Ø5x160,cannulated,inside Ø1.3,for cannulated screws Ø3.5	1
2.923.35	tap for quick coupling,Ø3.5x160,cannulated,inside Ø1.5	1
2.9436.25	hex screwdriver with handle for Ø2.7-4;cannulated;lØ1.3;wrench size 2.5;l 215mm	1
2.954.01	screw forceps,self holding	1
2.954.27A	holding sleeve for screws Ø2.7 to Ø3.5	1
2.955.12	cleaning wire Ø1.2 for cannulated instruments	1

Item no.	Item description	Quant.
Wires		
2.978.12	guide wire,Ø1.2x150,thread 5,st.st.	5

Implants - washers		
3.200.02	washer,outside Ø7 / inside Ø3,8,ti	5

Container		
19.185.00	perforated autoclavable container w.inset for instruments/cannul.screws Ø3.5 sd	1



Set no. 19.184.00  Internal hexagon