

Set Cannulated screws 3.5 st


self tapping; titanium; Set no. 19.181.00



Item no. 3.502.30



Item no. 3.602.30



Set no. 19.181.00  Internal hexagon

Item no.	Item description	Quant.
Implants – screws		
Partial thread - Ø3.5 mm; self tapping; titanium		
3.502.10	length 10 mm; thread 5	2
3.502.12	length 12 mm; thread 5	2
3.502.14	length 14 mm; thread 5	2
3.502.16	length 16 mm; thread 6	2
3.502.18	length 18 mm; thread 7	2
3.502.20	length 20 mm; thread 8	2
3.502.22	length 22 mm; thread 9	2
3.502.24	length 24 mm; thread 10	2
3.502.26	length 26 mm; thread 12	2
3.502.28	length 28 mm; thread 14	2
3.502.30	length 30 mm; thread 14	2
3.502.32	length 32 mm; thread 14	2
3.502.34	length 34 mm; thread 14	2
3.502.36	length 36 mm; thread 14	2
3.502.38	length 38 mm; thread 14	2
3.502.40	length 40 mm; thread 14	2
3.502.42	length 42 mm; thread 15	2
3.502.44	length 44 mm; thread 15	2
3.502.46	length 46 mm; thread 15	2
3.502.48	length 48 mm; thread 15	2
3.502.50	length 50 mm; thread 15	2
3.502.52	length 52 mm; thread 16	2
3.502.54	length 54 mm; thread 16	2
3.502.56	length 56 mm; thread 16	2

Item no.	Item description	Quant.
Implants – screws		
Full thread - Ø3.5 mm; self tapping; titanium		
3.602.10	length 10 mm	2
3.602.12	length 12 mm	2
3.602.14	length 14 mm	2
3.602.16	length 16 mm	2
3.602.18	length 18 mm	2
3.602.20	length 20 mm	2
3.602.22	length 22 mm	2
3.602.24	length 24 mm	2
3.602.26	length 26 mm	2
3.602.28	length 28 mm	2
3.602.30	length 30 mm	2
3.602.32	length 32 mm	2
3.602.34	length 34 mm	2
3.602.36	length 36 mm	2
3.602.38	length 38 mm	2
3.602.40	length 40 mm	2
3.602.42	length 42 mm	2
3.602.44	length 44 mm	2
3.602.46	length 46 mm	2
3.602.48	length 48 mm	2
3.602.50	length 50 mm	2
3.602.52	length 52 mm	2
3.602.54	length 54 mm	2
3.602.56	length 56 mm	2

Set Cannulated screws 3.5 st

self tapping; titanium; Set no. 19.181.00

Item no.	Item description	Quant.
3.502.58	length 58 mm; thread 16	2
3.502.60	length 60 mm; thread 16	2
3.502.62	length 62 mm; thread 17	2
3.502.64	length 64 mm; thread 17	2
3.502.66	length 66 mm; thread 17	2
3.502.68	length 68 mm; thread 17	2
3.502.70	length 70 mm; thread 17	2


Item no.	Item description	Quant.
3.602.58	length 58 mm	2
3.602.60	length 60 mm	2
3.602.62	length 62 mm	2
3.602.64	length 64 mm	2
3.602.66	length 66 mm	2
3.602.68	length 68 mm	2
3.602.70	length 70 mm	2

Instruments		
2.950.150	gauge for cannulated screws Ø 2.7 - 4.0 mm; length 150 mm	1
10.306.06	T-handle with quick coupling; length 93.5 mm; cannulated; inside Ø 2.5 mm	1
10.935.00	percutaneous drill bush system for cannulated screws Ø 3.5 mm	1
2.905.27	drill bit for quick coupling; Ø 2.7 x 160 mm; spiral length 30 mm; quadruple spiral; cannulated; inside Ø 1.35 mm	1
2.905.35	drill bit for quick coupling; Ø 3.5 x 160 mm; spiral length 30 mm; quadruple spiral; cannulated; inside Ø 1.35 mm	1
2.907.35	countersink; Ø 5.0 x 160 mm; cannulated; inside Ø 1.3 mm; for cannulated screws Ø 3.5 mm	1
2.923.35	tap for quick coupling; Ø 3.5 x 160 mm; cannulated; inside Ø 1.5 mm	1
2.9436.25	hex screwdriver with handle for Ø 2.7 - 4.0 mm; cannulated; inside Ø 1.3 mm; wrench size 2.5 mm; length 215 mm	1
2.954.01	screw forceps; self holding	1
2.954.27A	holding sleeve for screws Ø 2.7 - 3.5 mm	1
2.955.12	cleaning wire Ø 1.2 mm for cannulated instruments	
20.502.01	holding clamp; small	

Wires		
2.978.12	guide wire; Ø 1.2 x 150 mm; thread 5; stainless steel	5

Implants – washers		
3.200.02	washer; outside Ø 7.0 mm / inside Ø 3.8 mm; titanium	5

Container		
19.182.00	perforated autoclavable container with inset for instruments/cannulated screws Ø 3.5 mm	1



Set no. 19.181.00  Internal hexagon