

Set Cannulated screw 4.5 st sd


self tapping; self drilling; Titanium; Set no. 19.203.00



Item no. 3.434.30



Item no. 3.439.30



Set no. 19.203.00  Internal hexagon

Item no.	Item description	Quant.
Implants – screws		
Partial thread - Ø 4.5 mm; self drilling; back tapping; self tapping; titanium		
3.434.20	length 20 mm; thread 7	2
3.434.22	length 22 mm; thread 7	2
3.434.24	length 24 mm; thread 8	2
3.434.26	length 26 mm; thread 9	2
3.434.28	length 28 mm; thread 9	2
3.434.30	length 20 mm; thread 10	2
3.434.32	length 32 mm; thread 11	2
3.434.34	length 34 mm; thread 11	2
3.434.36	length 36 mm; thread 12	2
3.434.38	length 38 mm; thread 13	2
3.434.40	length 40 mm; thread 13	2
3.434.42	length 42 mm; thread 14	2
3.434.44	length 44 mm; thread 15	2
3.434.46	length 46 mm; thread 15	2
3.434.48	length 48 mm; thread 16	2
3.434.50	length 50 mm; thread 16	2
3.434.52	length 52 mm; thread 17	2
3.434.54	length 54 mm; thread 19	2
3.434.56	length 56 mm; thread 19	2
3.434.58	length 58 mm; thread 20	2
3.434.60	length 60 mm; thread 20	2
3.434.62	length 62 mm; thread 21	2
3.434.64	length 64 mm; thread 21	2
3.434.66	length 66 mm; thread 22	2
3.434.68	length 68 mm; thread 23	2
3.434.72	length 72 mm; thread 24	2

Item no.	Item description	Quant.
Implants – screws		
Full thread - Ø 4.5 mm; self drilling; self tapping; titanium		
3.439.20	length 20 mm	2
3.439.22	length 22 mm	2
3.439.24	length 24 mm	2
3.439.26	length 26 mm	2
3.439.28	length 28 mm	2
3.439.30	length 30 mm	2
3.439.32	length 32 mm	2
3.439.34	length 34 mm	2
3.439.36	length 36 mm	2
3.439.38	length 38 mm	2
3.439.40	length 40 mm	2
3.439.42	length 42 mm	2
3.439.44	length 44 mm	2
3.439.46	length 46mm	2
3.439.48	length 48 mm	2
3.439.50	length 50 mm	2
3.439.52	length 52 mm	2
3.439.54	length 54 mm	2
3.439.56	length 56 mm	2
3.439.58	length 58 mm	2
3.439.60	length 60 mm	2
3.439.62	length 62 mm	2
3.439.64	length 64 mm	2
3.439.66	length 66 mm	2
3.439.68	length 68 mm	2
3.439.72	length 72 mm	2

Set Cannulated screw 4.5 st sd


self tapping; self drilling; Titanium; Set no. 19.203.00

Item no.	Item description	Quant.
Instruments		
2.951.150	gauge for cannulated screws Ø4.5,length 150	1
10.306.06	T-handle with quick coupling,length 93.5,cannulated,inside Ø2.5	1
10.945.00	percutaneous drill bush system for cannulated screws Ø4.5	1
20.502.01	holding clamp,small	1
2.905.32	drill bit f.quick coupling,Ø3.2x160,spiral ln 60,quadr.spiral,cann.,inside Ø1.75	1
2.905.45	drill bit f.quick coupling,Ø4.5x160,spiral ln 60,quadr.spiral,cann.,inside Ø1.8	1
2.907.45	countersink,Ø6.5x160,cannulated,inside Ø1.75,for cannulated screws Ø4.5	1
2.923.45	tap with T-handle,Ø4.5x135,cannulated,inside Ø1.75	1
2.9436.35	hex screwdriver with handle for Ø4.5-7;cannulated;lØ2.2;wrench size 3.5;l 215mm	1
2.954.01	screw forceps,self holding	1
2.954.45A	holding sleeve for screws Ø4.5	1
2.955.16	cleaning wire Ø1.6 for cannulated instruments	1

Item no.	Item description	Quant.
Wires		
2.978.16	guide wire,Ø1.6x150,thread 5,st.st.	5

Implants - washers		
3.200.40	washer,outside Ø10 / inside Ø5,ti	5

Container		
19.205.00	perforated autoclavable container with inset for instruments/cannul.screws Ø4.5	1



Set no. 19.203.00  Internal hexagon