

Set Cannulated screws 6.5 st sd IH T


self tapping; self drilling; Titanium; TiOB®; Set no. 19.210.00



Item no. 3.520.70



Item no. 3.523.70



Set no. 19.210.00 Internal hexagon

Item no.	Item description	Quant.
Implants – screws		
Cannulated screw; partial thread 16 - Ø 6.5 mm; self drilling; back tapping; self tapping; titanium; TiOB®		
3.520.30	length 30 mm	2
3.520.35	length 35 mm	2
3.520.40	length 40 mm	2
3.520.45	length 45 mm	2
3.520.50	length 50 mm	2
3.520.55	length 55 mm	2
3.520.60	length 60 mm	2
3.520.65	length 65 mm	2
3.520.70	length 70 mm	2
3.520.75	length 75 mm	2
3.520.80	length 80 mm	2
3.520.85	length 85 mm	2
3.520.90	length 90 mm	2
3.520.95	length 95 mm	2
3.520.100	length 100 mm	2
3.520.105	length 105 mm	2
3.520.110	length 110 mm	2
3.520.115	length 115 mm	2
3.520.120	length 120 mm	2
3.520.125	length 125 mm	2
3.520.130	length 130 mm	2

Item no.	Item description	Quant.
Implants – screws		
Cannulated screw; partial thread 32 - Ø 6.5 mm; self drilling; back tapping; self tapping; titanium; TiOB®		
3.523.45	length 45 mm	2
3.523.50	length 50 mm	2
3.523.55	length 55 mm	2
3.523.60	length 60 mm	2
3.523.65	length 65 mm	2
3.523.70	length 70 mm	2
3.523.75	length 75 mm	2
3.523.80	length 80 mm	2
3.523.85	length 85 mm	2
3.523.90	length 90 mm	2
3.523.95	length 95 mm	2
3.523.100	length 100 mm	2
3.523.105	length 105 mm	2
3.523.110	length 110 mm	2
3.523.115	length 115 mm	2
3.523.120	length 120 mm	2
3.523.125	length 125 mm	2
3.523.130	length 130 mm	2

Information

TiOB® is a surface treatment that brings optimized properties.



Königsee Implantate GmbH
OT Aschau | Am Sand 4
07426 Allendorf | Germany

Fon +49 36738 498-560
Fax +49 36738 498-569

Mail export@koenigsee-implantate.de

For further information please visit:
www.koenigsee-implantate.com

HeBo | Edition 05 | 2022-12
All pictured products are CE compliant.

Set Cannulated screws 6.5 st sd IH T


self tapping; self drilling; Titanium; TiOB®; Set no. 19.210.00

Item no.	Item description	Quant.
Instruments		
2.952.230	gauge for cannulated screws Ø 6,5 mm and Ø 7,0 mm, length 230 mm	1
10.306.06	T-handle with quick coupling, length 93.5 mm, cannulated, inside Ø 2.5 mm	1
10.965.00	percutaneous drill bush system for cannulated screws Ø 6.5 mm	1
2.907.70	countersink, Ø8x160mm, cutter ln 19.5 mm, cannulated, inside Ø 2.1 mm, for cannulated screws Ø 6.5 mm	1
2.908.40	drill bit w.triangular shank, Ø 4x230 mm, spiral ln 40 mm, cannulated, inside Ø 2.1 mm, triple sp.	1
2.942.02	hex screwdriver with T-handle for screw Ø 4.5 mm until Ø 7.0 mm, cannulated, inside Ø 2.2 mm	1
2.954.01	screw forceps, self holding	1
2.954.65A	holding sleeve for screws Ø 6.5 mm and Ø 7,0 mm	1
2.955.20	cleaning wire Ø 2.0 mm for cannulated instruments	1
2.960.03	parallel targeting device for guide wire Ø 2.0 mm	1

Item no.	Item description	Quant.
Wires		
2.979.20	guide wire, Ø 2x230 mm, thread 5, stainless steel	5

Implants - washers		
3.200.04	washer, outside Ø 16 mm / inside Ø 6.7 mm, titan	5

Container		
19.211.00	perforated autoclavable container with inset for instruments/cannulated screws Ø 6.5 mm	1



Set no. 19.210.00  Internal hexagon