


# Set FRAKO 3.0/4.0 DrillFix 2.0 INSTRUMENTS

Set no. 19.335.001



Set no. 19.335.001  Internal hexagon

Item no.	Item description	Quant.
<b>Instruments</b>		
2.906.30	Pilot Drill Bit Frako qc 2.1; DL 8 mm; Burr. 3.0; cann. 1.3; length 60 mm	1
2.906.31	Pilot Drill Bit Frako qc 2.1; DL 18 mm; Burr. 3.0; cann. 1.3; length 60 mm	1
10.3066.01	handle with quick coupling, length 146 mm, cannulated, inside Ø 2.5 mm; PPSU	1
2.950.01	gauge for fracture compression screw; length 150 mm; MR 40 mm	1
2.9436.20	hex screwdriver with handle for fracture compression; cann.; IØ 1.3 mm; wrench size 2; length 215 mm	1
2.9406.127	Screwdriver with handle drill fix screw	1
2.956.42	tissue protecting drill sleeve, Ø 6.0 x 30 mm, inside Ø 4.2 mm	1

Item no.	Item description	Quant.
<b>Wires</b>		
2.982.12	guide wire, Ø 1.2 x 150 mm, thread 10, stainless steel	5
6.031.12	Kirschner wire with trocar point and round end, Ø 1.2 x 150 mm, stainless steel	3

<b>Container</b>		
19.336.00	perforated autoclavable container with inset for instruments and fracture compression screw	1

#### Can be combined with:

19.335.00	FRAKO 3.0/4.0 DrillFix 2.0 SET
19.335.00PS+	FRAKO 3.0/4.0 SCREWS sterile