



## All advantages at a glance

- Self-drilling and self-tapping screws
- **Cannulated** fracture compression screws enable a **minimally invasive** surgical technique
- Drill-Fix forefoot screw designed as **Break-off** screw
- **Storage space saving**, because the screws are available in one storage system



## Contact details

If you are interested in the  
**Fracture compression screw / Drill-fix forefoot screw**  
please feel free to contact us or your local sales  
representatives.



### Königsee Implantate GmbH

OT Aschau | Am Sand 4  
07426 Allendorf  
Germany

[export@koenigsee-implantate.de](mailto:export@koenigsee-implantate.de)

Fon: + 49 36738 498 - 560


Fax: + 49 36738 498 - 89 560

All implants are available on request as non-sterile  
packed or sterile packed.

[www.koenigsee-implantate.de](http://www.koenigsee-implantate.de)

## FRAKO & DRILL-FIX

Forefoot treatment / Fracture compression



## Available screw lengths

### Drill-fix forefoot screws


Ø 2.0 mm, titanium

thread length 11 mm; total length 26 mm


thread length 12 mm; total length 27 mm

thread length 13 mm; total length 28 mm

thread length 14 mm; total length 29 mm



Fracture compression screw



Drill-fix forefoot screw

## Available screw lengths

### Fracture compression screws

Ø 3.0 mm / Ø 4.0 mm,

self drilling; self tapping;

cannulated; inside Ø 1.3 mm; titanium

total length 10 mm; thread length 4 mm

total length 12 mm; thread length 4 mm

total length 14 mm; thread length 4 mm

total length 16 mm; thread length 6 mm

total length 18 mm; thread length 6 mm

total length 20 mm; thread length 6 mm

total length 22 mm; thread length 6 mm


total length 24 mm; thread length 6 mm

total length 26 mm; thread length 6 mm

total length 28 mm; thread length 6 mm

total length 30 mm; thread length 6 mm

## Set FRAKO and DRILL-FIX



Set fracture compression screws and drill-fix forefoot screws



Self holding screw forceps



Cannulated gauge



Hex screwdriver

Tissue protecting drill sleeve

