

Set FRAKO 3.0/4.0 DrillFix 2.0


Set no. 19.335.00



Item no. 3.123.20



Item no. 3.127.12



Set no. 19.335.00 Internal hexagon

Item no.	Item description	Quant.
Implants		
Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm; self drilling; self tapping; Inside Ø 1.3 mm; titanium		
3.123.10	length 10 mm; thread 4	3
3.123.12	length 12 mm; thread 4	3
3.123.14	length 14 mm; thread 4	3
3.123.16	length 16 mm; thread 6	3
3.123.18	length 18 mm; thread 6	3
3.123.20	length 20 mm; thread 6	3
3.123.22	length 22 mm; thread 6	3
3.123.24	length 24 mm; thread 6	3
3.123.26	length 26 mm; thread 6	3
3.123.28	length 28 mm; thread 6	3
3.123.30	length 30 mm; thread 6	3
Drill Fix screw; titanium		
3.127.11	length 11 mm; total length 26 mm	3
3.127.12	length 12 mm; total length 27 mm	3
3.127.13	length 13 mm; total length 28 mm	3
3.127.14	length 14 mm; total length 29 mm	3

Instruments		
2.906.30	Pilot Drill Bit Frako qc 2.1; DL 8; Burr. 3.0; cann. 1.3; length 60 mm	1
2.906.31	Pilot Drill Bit Frako qc 2.1; DL 18; Burr. 3.0; cann. 1.3; length 60 mm	1
10.3066.01	handle with quick coupling, length 146 mm, cannulated, inside Ø 2.5 mm; PPSU	1

Instruments		
2.950.01	gauge for fracture compression screw; length 150 mm; measuring range 40 mm	1
2.9436.20	hex screwdriver with handle for fracture compression; cann.; lØ 1.3 mm; wrench size 2; length 15 mm	1
2.9406.127	Screwdriver with handle drill fix screw	1
2.956.42	tissue protecting drill sleeve, Ø 6.0 x 30 mm, inside Ø 4.2 mm	1
2.954.01	screw forceps, self holding	1


Wires		
2.982.12	guide wire, Ø 1.2 x 150 mm, thread 10, stainless steel	5
6.031.12	Kirschner wire with trocar point and round end, Ø 1.2 x 150 mm, stainless steel	3

Container		
19.336.00	perforated autoclavable container with inset for instruments and fracture compression screw	1

All implants are available on request as non-sterile packed or sterile packed.

Set FRAKO 3.0/4.0 DrillFix 2.0

Set no. 19.335.00



Set no. 19.335.00  Internal hexagon