

# Set FRAKO 3.0/4.0 DrillFix 2.0 / SCREWS sterile

for forefoot treatment; titanium


non-sterile packaging Set no. 19.335.00 / sterile packaging Set no. 19.335.00PS+



Item no. 3.123.20  
Item no. S3.123.20




Item no. 3.127.12  
Item no. S3.127.12



All listed implants are available individually in non-sterile and sterile packaging.

Implants in sets, supplied **non-sterile** Set no. 19.335.00

Implants in sets, supplied **sterile** Set no. 19.335.00PS+

 **Internal hexagon**  
Set no. 19.335.00  
non-sterile packaging

Item no.	Item no. sterile	Item description	Quant.
<b>Implants – screws</b>			
<b>Fracture compression screw 3.0/4.0 st sd Cann.</b>			
3.123.10	S3.123.10	Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 10 mm; thread length 4.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium	3
3.123.12	S3.123.12	Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 12 mm; thread length 4.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium	3
3.123.14	S3.123.14	Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 14 mm; thread length 4.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium	3
3.123.16	S3.123.16	Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 16 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium	3
3.123.18	S3.123.18	Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 18 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium	3
3.123.20	S3.123.20	Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 20 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium	3
3.123.22	S3.123.22	Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 22 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium	3
3.123.24	S3.123.24	Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 24 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium	3
3.123.26	S3.123.26	Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 26 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium	3
3.123.28	S3.123.28	Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 28 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium	3
3.123.30	S3.123.30	Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 30 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium	3

<b>Drill Fix screw</b>			
3.127.11	S3.127.11	DrillFix screw; length 11 mm; total length 26 mm; titanium	3
3.127.12	S3.127.12	DrillFix screw; length 12 mm; total length 27 mm; titanium	3
3.127.13	S3.127.13	DrillFix screw; length 13 mm; total length 28 mm; titanium	3
3.127.14	S3.127.14	DrillFix screw; length 14 mm; total length 29 mm; titanium	3

# Set FRAKO 3.0/4.0 DrillFix 2.0 / SCREWS sterile

for forefoot treatment; titanium

non-sterile packaging Set no. 19.335.00 / sterile packaging Set no. 19.335.00PS+

Item no.	Item no. sterile	Item description	Quant.
<b>Instruments</b>			
2.906.30	S2.906.30	Pilot drill; fracture compression screw; Ø 2.1 mm; drill length 8.0 mm; burr Ø 3.0 mm; drilled inner Ø 1.3 mm; length 60 mm; quick coupling	1
2.906.31	S2.906.31	Pilot drill; fracture compression screw; Ø 2.1 mm; drill length 18 mm; burr Ø 3.0 mm; drilled inner Ø 1.3 mm; length 60 mm; quick coupling	1
10.3066.01	—	Handle with quick coupling; perforated; inner Ø 2.5 mm; length 146 mm; PPSU	1
2.950.01	—	Gauge for fracture compression screw; cannulated; measuring range 40 mm; length 150 mm	1
2.9436.20	—	Hexagon screwdriver; for hexagon socket 2.0; for fracture compression screw; perforated; inner Ø 1.3 mm; length 215 mm	1
2.9406.127	—	Screwdriver; for DrillFix screws; length 215 mm; PPSU handle	1
2.956.42	—	Drill and tissue protection sleeve; two-piece; Ø 6.0 mm x 30 mm; inner Ø 4.2 mm	1
2.954.01	—	Tweezer for screws self-holding	1


<b>Wires</b>			
2.982.12	—	Guide wire; Ø 1.2 mm x 150 mm; thread 10 mm; stainless steel	5
6.031.12	S6.031.12	Kirschner wire with trocar tip and round end; Ø 1.2 mm x 150 mm; stainless steel	3


<b>Container</b>			
19.336.00	—	Perforated container with insert for instruments; fracture compression screws and DrillFix screws	1

# Set FRAKO 3.0/4.0 DrillFix 2.0 / SCREWS sterile

for forefoot treatment; titanium

non-sterile packaging Set no. 19.335.00 / sterile packaging Set no. 19.335.00PS+



 Internal hexagon  
Set no. 19.335.00  
non-sterile packaging