


Fracture compression screw 3.0/4.0

for forefoot treatment / Dens screwing; titanium

non-sterile packaging / sterile packaging



Item no. 3.123.20















Item no. S3.123.20

| Item no. | EAN | Item description | Item no. sterile | EAN sterile |
|--|-----|--|------------------|-------------|
| Fracture compression screw 3.0/4.0 st sd Cann. – for forefoot treatment | | | | |
| 3.123.10 | | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 10 mm; thread length 4.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.10 | |
| 3.123.12 | | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 12 mm; thread length 4.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.12 | |
| 3.123.14 | | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 14 mm; thread length 4.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.14 | |
| 3.123.16 | | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 16 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.16 | |
| 3.123.18 | | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 18 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.18 | |
| 3.123.20 | | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 20 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.20 | |
| 3.123.22 | | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 22 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.22 | |
| 3.123.24 | | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 24 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.24 | |
| 3.123.26 | | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 26 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.26 | |
| 3.123.28 | | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 28 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.28 | |
| 3.123.30 | | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 30 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.30 | |

Fracture compression screw 3.0/4.0

for forefoot treatment / Dens screwing; titanium

non-sterile packaging / sterile packaging

| Item no. | EAN | Item description | Item no. sterile | EAN sterile |
|--|--|--|------------------|--|
| Fracture compression screw 3.0/4.0 st sd Cann. – for Dens screwing | | | | |
| 3.123.32 |  4 047075 088805 | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 32 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.32 |  4 047075 191536 |
| 3.123.34 |  4 047075 088869 | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 34 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.34 |  4 047075 191543 |
| 3.123.36 |  4 047075 088876 | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 36 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.36 |  4 047075 191550 |
| 3.123.38 |  4 047075 088883 | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 38 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.38 |  4 047075 191557 |
| 3.123.40 |  4 047075 088890 | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 40 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.40 |  4 047075 191564 |
| 3.123.42 |  4 047075 088706 | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 42 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.42 |  4 047075 191581 |
| 3.123.44 |  4 047075 088560 | Fracture compression screw; Ø 3.0 mm/Ø 4.0 mm x 44 mm; thread length 6.0 mm; self tapping; self drilling; internal Ø 1.3 mm; internal hexagon 2.0 mm; titanium | S3.123.44 |  4 047075 191598 |