


# Set Instruments TIS screw 7.5

Set no. 19.075.001



Set no. 19.075.001  Internal hexagon

Item no.	Item description	Quant.
Instruments		
2.908.50	drill bit with coupling; Ø 5 x 300 mm; spiral length 55 mm; cannulated; Ø 2.9 mm; four-spiral	1
2.9435.75	hex screwdriver with handle for Ø 7.5 mm; cannulated; inside Ø 2.9 mm; wrench size 4; length 324 mm	1
2.952.200	gauge for cannulated screws Ø 7.5 mm, measuring range 200 mm, length 275 mm	1
20.502.01	holding clamp; small	1

Item no.	Item description	Quant.
Implants – washers		
3.200.04	washer; outside Ø 16 mm / inside Ø 6.7 mm; titanium	5
Wires		
2.980.28	guide wire; Ø 2.8 x 300 mm; thread 5; stainless steel	5
2.980.28/450	guide wire; Ø 2.8 x 450 mm; thread 5; stainless steel	5
Container		
19.075.01	Perforated autoclavable container for instruments TIS	1

Can be combined with:

19.681.00	TIS screws 7.5 st IH T SET
19.681.00/1	TIS screws 7.5 st IH T upgrade SCREWS
19.681.00PS	TIS screws 7.5 st IH T SCREWS sterile MIN
19.681.00/1PS+	TIS screws 7.5 st IH T upgrade SCREWS sterile
19.681.00PS+	TIS screws 7.5 st IH T SCREWS sterile