


Set Instruments LF IH

Large Fragment; for screws Ø 4.5 / 5.5 / 6.0 / 6.5 mm; Set no. 19.045.001



Set no. 19.045.001  Internal hexagon


Item no.	Item description	Quant.
Instruments		
2.904.09	Drill bit qc 3.2; Spiral 60 mm; 2-flute; length 195 mm	1
2.904.23	Drill bit qc 4.0; Spiral 60 mm; 2-flute; length 220 mm	1
2.9406.35	hex screwdriver with handle for screws Ø 4.5 mm - Ø 7.0 mm; AF 3.5; length 215 mm	1
2.953.90	gauge with clasp for screws Ø 4.5 mm; Ø 6.0 mm; measuring range 90 mm	1
2.977.02	drill guide for angle stable screwing; length 70 mm; large fragment	2
2.977.15	drill guide for angle stable screwing 5.5; length 70 mm	2
2.9771.45	drill guide for angle stable compression hole; large fragment	1
10.502.32	universal drill guide; Ø 3.2 mm	1

Item no.	Item description	Quant.
Wires		
6.031.16	Kirschner wire with trocar point and round end; Ø 1.6 x 150 mm; stainless steel	5
6.031.18	Kirschner wire with trocar point and round end; Ø 1.8 x 150 mm; stainless steel	5

Container		
19.045.01	Perforated autoclavable container for instruments large fragment; LF	1

Set Instruments LF IH

Large Fragment; for screws Ø 4.5 / 5.5 / 6.0 / 6.5 mm; Set no. 19.045.001



Can be combined with:

19.103.00	Prox.Humerus plate 6.0/4.5 as IH S SET
19.104.00	Prox.Humerus plate 6.0/4.5 as IH T SET
19.600.00	Distal Femur plate 4.5/5.5 as IH S SET
19.600.00PS+	Distal Femur plate 4.5/5.5 as IH S IMPLANTS sterile
19.602.00	Distal Femur plate 4.5/5.5 as IH T SET
19.602.00PS+	Distal Femur plate 4.5/5.5 as IH T IMPLANTS sterile
19.404.00	T Tibia head plate 4.5/6.0 as IH S SET
19.404.00PS+	T Tibia head plate 4.5/6.0 as IH S IMPLANTS sterile
19.408.00	T Tibia head plate 4.5/6.0 as IH T SET
19.408.00PS+	T Tibia head plate 4.5/6.0 as IH T IMPLANTS sterile
19.402.00	L Tibia head plate 4.5/6.0 as IH S SET
19.402.00PS+	L Tibia head plate 4.5/6.0 as IH S IMPLANTS sterile
19.406.00	L Tibia head plate 4.5/6.0 as IH T SET
19.406.00PS+	L Tibia head plate 4.5/6.0 as IH T IMPLANTS sterile
19.407.00	T & L Tibia head plate 4.5/6.0 as IH T SET
19.409.00	Med. Tibia head plate 4.5/6.0 as IH T SET
19.409.00PS+	Med. Tibia head plate 4.5/6.0 as IH T IMPLANTS sterile
19.516.00	CPL plate 4.5 as IH S SET
19.161.00	LF plate S PLATES
19.162.00	LF plate T PLATES
19.512.00	Low Contact plate 4.5/6.0 as IH S SET
19.514.00	Low Contact plate 4.5/6.0 as IH T SET

Set no. 19.045.001  Internal hexagon