


Set Instruments 2.4 SL

for screws Ø 2.4 mm; Set no. 19.024.001



Set no. 19.024.001  Six lobe

Item no.	Item description	Quant.
Instruments		
2.904.33	drill bit for quick coupling; Ø1.8 x 105 mm; spiral length 30 mm; double spiral	1
2.953.24	gauge with clasp for screws Ø 2.4 mm; measuring range 60 mm	1
2.9775.24	drill guide; Ø 2.4 mm double cone; angle-stable	2
2.93886.08	screwdriver hexalobe size 8; PPSU-handle; length 180 mm	1
10.500.24	drill guide with handle for standard screwing 2.4; length 45 mm	1

Item no.	Item description	Quant.
Wires		
6.031.15	Kirschner wire with trocar point and round end; Ø 1.5 x 150 mm; stainless steel	5
Container		
19.024.01	Perforated autoclavable container for Instruments 2.4	1

Can be combined with: 19.117.00 Dist. radius plate 2.4 vas SL T SET
 19.117.00PS Dist. radius plate 2.4 vas SL T IMPLANTS sterile MIN
 19.117.00PS+ Dist. radius plate 2.4 vas SL T IMPLANTS sterile