


Set Instruments 2.7 IH

for screws Ø 2.7 mm; Set no. 19.027.001



Set no. 19.027.001  Internal hexagon

Item no.	Item description	Quant.
Instruments		
2.904.03	drill bit for quick coupling; Ø 2x105 mm; spiral length 30 mm; double spiral	1
2.9406.25	hex screwdriver with handle for screws Ø 2.7 mm - Ø 4.0 mm; AF 2.5; conical	1
2.953.65	gauge with clasp for screws with conical head thread Ø 2.7; measuring range 60	1
2.977.19	drill guide w.handle for variable angle stable screwing 2.7; length 42	1
2.977.11/3	drill guide for angle stable screwing 2.7; length 44 mm; single thread with knurl	2

Item no.	Item description	Quant.
Wires		
6.031.16	Kirschner wire with trocar point and round end; Ø 1.6 x 150 mm; stainless steel	5
Container		
19.027.01	Perforated autoclavable container for instruments 2.7	1

Optionally available:

- 10.180.05 – guide block for distal radius plate 5.9823.xx; right
- 10.181.05 – guide block for distal radius plate 5.9833.xx; left
- 10.180.07 – guide block for distal radius plate 5.9843.xx; right
- 10.181.07 – guide block for distal radius plate 5.9853.xx; left

Can be combined with:

19.116.00	Dist. radius plate 2.7 as IH T SET
19.116.00PS	Dist. radius plate 2.7 vas IH T SET
19.116.00PS+	Dist. radius plate 2.7 as IH T IMPLANTS sterile MIN
19.1163.00	Dist. radius plate 2.7 vas IH T IMPLANTS sterile MIN
19.1163.00PS	Dist. radius plate 2.7 as IH T IMPLANTS sterile
19.1163.00PS+	Dist. radius plate 2.7 vas IH T IMPLANTS sterile
5.3603.xx	Clavicle plate with conical thread 2,7; R; vas; T
5.3613.xx	Clavicle plate with conical thread 2,7; L; vas; T