


Set Instruments 2.7 SL

for screws Ø 2.7 mm; Set no. 19.027.10I



Set no. 19.027.10I Six lobe

Item no.	Item description	Quant.
Instruments		
2.904.03	drill bit for quick coupling; Ø 2.0 x 105 mm; spiral length 30 mm; double spiral	1
2.9386.10	screwdriver with handle for six-lobe 10	1
2.953.65	gauge with clasp for screws Ø 2.7 mm; measuring range 60 mm	1
2.977.19	drill guide with handle for variable angle stable screwing 2.7; length 42 mm	1
2.977.11/3	drill guide for angle stable screwing 2.7; length 44 mm; single thread with knurl	2

Item no.	Item description	Quant.
Wires		
6.031.16	Kirschner wire with trocar point and round end; Ø 1.6 x 150 mm; stainless steel	5
Container		
19.027.01	Perforated autoclavable container for instruments 2.7	1

Optionally available:

10.180.05 – guide block for distal radius plate 5.9823.xx; right
 10.181.05 – guide block for distal radius plate 5.9833.xx; left
 10.180.07 – guide block for distal radius plate 5.9843.xx; right
 10.181.07 – guide block for distal radius plate 5.9853.xx; left

Can be combined with: 19.1161.00 Dist. radius plate 2.7 vas SL T SET
 19.1161.00PS Dist. radius plate 2.7 vas SL T IMPLANTS sterile MIN
 19.1161.00PS+ Dist. radius plate 2.7 vas SL T IMPLANTS sterile