


Set Instruments SF IH

for screws Ø 3.5/4.0; Set no. 19.035.001



Set no. 19.035.001 - Level 1



Set no. 19.035.001 - Level 2



Set no. 19.035.001  Internal hexagon

LEVEL 1

Item no.	Item description	Quant.
Instruments		
2.9406.25	hex screwdriver with handle for screws Ø 2.7 – Ø 4.0; AF2.5; conical	1
2.904.05	drill bit for quick coupling; Ø 2.5 x 105; spiral length 50; double spiral	1
2.904.22	drill bit for quick coupling; Ø 2.5 x 155; spiral length 60; double spiral	1
10.500.27A	drill guide with handle for small fragment; handle length 55	1
2.904.25	drill bit for quick coupling; Ø 3.0 x 155; spiral length 60; double spiral	1
2.977.16	drill guide for angle stable screwing 4.0; length 68	2
2.977.17	drill guide w.handle for variable angle stable screwing 4.0; length 68	1

Wires

6.031.16	Kirschner wire with trocar point and round end; Ø 1.6 x 150; st.st.	5
----------	---	---

Optionally available – Level 1:

- 2.934.25 – torque screwdriver handle; 2.2 newtonmeter
- 2.941.25A – hex screwdriver inset f.quick coupling f.screw Ø 2.7 – Ø 4.0; wrench width 2.5; conic
- 2.977.12 – drill guide with handle for variable angle stable screwing 3.5 with head Ø 5.0
- 2.977.13 – drill guide for angle stable screws 3.5 with head Ø 5.0; length 44
- 2.904.10 – drill bit for quick coupling; Ø3.5x110; spiral length 40; double spiral
- 10.520.30 – double drill guide Ø2.5/Ø3.5; spring bushing

Optionally available – Level 2:

- 2.9771.35 – drill guide for angle stable compression hole; small fragment
- 2.911.06 – bending lever for reco plates small fragment

LEVEL 2

Item no.	Item description	Quant.
Instruments – Level 2; part 1		
2.977.04	drill guide for angle stable screwing; length 68; small fragment	2
2.977.09	drill guide w.handle for variable angle stable screwing; small fragment; length 68	1
2.953.70	gauge w. clasp for screws Ø 3.5; Ø 4.0; measuring range 70	1
10.179.10	guide block for proximal small fragment humerus plates; 8 holes in head; flat	1
10.173.00	guide block for proximal small fragment humerus plates w. 5 headholes	1
10.179.15	aiming block f. prox.hum. pl.5.07913.8X	1
10.173.01	sleeve for guide block for proximal small fragment humerus plates	2
10.173.05	screw driver for guide block	1

Instruments – Level 2; part 2

2.977.01	drill guide for angle stable screwing; length 45; small fragment	2
2.977.08	drill guide w.handle for variable angle stable screwing; small fragment; length 45	1
2.953.60	gauge w. clasp for screws Ø 3.5; Ø 4.0; measuring range 60	1

Container

19.530.00	perforated autoclavable container with inset f.instruments-modul small fragment	1
-----------	---	---

Set Instruments SF IH


for screws Ø 3.5/4.0; Set no. 19.035.001

Can be combined with:

19.478.00	Clavicle plate 3.5 as IH T SET	19.367.00	Prox.Humerus plate FIVE vas IH T SET
19.478.10PS	Clavicle plate lateral as IH T IMPLANTS sterile MIN	19.367.00PS	Prox.Humerus plate FIVE vas IH T IMPLANTS sterile MIN
19.478.00PS	Clavicle plate 3.5 as IH T IMPLANTS sterile MIN	19.367.00PS+	Prox.Humerus plate FIVE vas IH T IMPLANTS sterile
19.478.20PS	Clavicle shaft plate as IH T IMPLANTS sterile MIN	19.3673.00	Prox.Humerus plate EIGHT vas IH T SET
19.478.00PS+	Clavicle plate 3.5 as IH T IMPLANTS sterile	19.3673.00PS	Prox.Humerus plate EIGHT vas IH T IMPLANTS sterile MIN
19.478.30PS	Clavicle hook plate as IH T IMPLANTS sterile MIN	19.3673.00PS+	Prox.Humerus plate EIGHT vas IH T IMPLANTS sterile
19.478.40PS	Clavicle plate anterior as IH T IMPLANTS sterile MIN	19.3683.00	Prox.Humerus plate EIGHT vas Corti IH T SET
19.479.00	Clavicle plate 3.5 vas IH T SET	19.3683.00PS	Prox.Humerus plate EIGHT vas Corti IH T IMPLANTS sterile MIN
19.479.10PS	Clavicle plate lateral vas IH T IMPLANTS sterile MIN	19.3683.00PS+	Prox.Humerus plate EIGHT vas Corti IH T IMPLANTS sterile
19.479.00PS	Clavicle plate 3.5 vas IH T IMPLANTS sterile MIN	19.3751.00	VaWiKo® Prox.humerus plate 3.5 IH SET
19.479.20PS	Clavicle shaft plate vas IH T IMPLANTS sterile MIN	19.3751.00PS	VaWiKo® Prox.humerus 4.0/3.5 IH IMPLANTS sterile MIN
19.479.00PS+	Clavicle plate 3.5 vas IH T IMPLANTS sterile	19.3751.00PS+	VaWiKo® Prox.humerus plate 4.0/3.5 IH IMPLANTS sterile
19.479.30PS	Clavicle hook plate vas IH T IMPLANTS sterile MIN	19.3752.00	VaWiKo® Prox.humerus 4.0/3.5 IH Corti SET
19.364.00	Prox.Humerus plate FIVE as IH T SET	19.3752.00PS	VaWiKo® Prox.humerus 4.0/3.5 IH Corti IMPLANTS sterile MIN
19.364.00PS	Prox.Humerus plate FIVE as IH T IMPLANTS sterile MIN	19.3752.00PS+	VaWiKo® Prox.humerus 4.0/3.5 IH Corti IMPLANTS sterile
19.364.00PS+	Prox.Humerus plate FIVE as IH T IMPLANTS sterile	19.3753.00	VaWiKo® Prox.humerus 4.0/3.5 IH Cancell SET
		19.3753.00PS	VaWiKo® Prox.humerus 4.0/3.5 IH IMPLANTS sterile MIN
		19.3753.00PS+	VaWiKo® Prox.humerus plate 4.0/3.5 IH IMPLANTS sterile



Set no. 19.035.001 - Level 1



Set no. 19.035.001 - Level 2