


# Set Instruments SF 3.5 IH Head 5.0

for screws Ø 3.5; Head 5.0; Set no. 19.035.10I



Set no. 19.035.10I  Internal hexagon

Item no.	Item description	Quant.
Instruments		
2.904.06	Drill bit qc2.5; Spiral 60 mm; 2-flute; length 125 mm; blue	1
2.977.12	drill guide with handle for variable angle stable screwing 3.5 with head Ø 5.0 mm	1
2.977.13	drill guide for angle stable screws 3.5 with head Ø 5.0 mm; length 44 mm	2
2.9406.25	hex screwdriver with handle for screws Ø 2.7 mm - Ø 4.0 mm; AF 2.5; conical	1
2.953.60	gauge with clasp for screws Ø 3.5 mm; Ø 4.0 mm; measuring range 60 mm	1

Item no.	Item description	Quant.
Wires		
6.031.10	Kirschner wire with trocar point and round end; Ø 1,0 x 150 mm; stainless steel	5

Container		
19.035.11	Perforated autoclavable container Instruments 3.5 Head Ø 5.0 mm IH	1

Can be combined with:	19.671.010	Patella plate 3.5 vas IH T SET
	19.671.010PS+	Patella plate 3.5 vas IH T IMPLANTS sterile
	19.524.00	Hintfoot plate 3.5 vas IH T SET
	19.524.00PS+	Hintfoot plate 3.5 vas IH T IMPLANTS sterile