

# Angle plates

## stainless steel



Item no. 4.676.01



Item no. 4.636.04











Item no. 4.656.022

| Item no.   | Item description   | EAN             |
|--|--|-----------------|
| <b>Adolescent plate intertroch. 90° 4.5 – plate thickness 4.0 mm</b> |  |                 |
| 4.676.01   | intertroch.plate for adolescents 90° 4.5 mm; 3 holes; arc depth 10 mm; blade 40 mm; stainless steel                  | 4 047075 053445 |
| 4.676.02   | intertroch.plate for adolescents 90° 4.5 mm; 3 holes; arc depth 10 mm; blade 50 mm; stainless steel                  | 4 047075 053452 |
| 4.676.03   | intertroch.plate for adolescents 90° 4.5 mm; 3 holes; arc depth 15 mm; blade 40 mm; stainless steel                  | 4 047075 053513 |
| 4.676.04   | intertroch.plate for adolescents 90° 4.5 mm; 3 holes; arc depth 15 mm; blade 50 mm; stainless steel                  | 4 047075 053520 |
| 4.676.05   | intertroch.plate for adolescents 90° 4.5 mm; 3 holes; arc depth 15 mm; blade 60 mm; stainless steel                  | 4 047075 053544 |
| 4.676.23   | intertroch.plate for adolescents 90° 4.5 mm; 3 holes; arc depth 10 mm; blade 55 mm; stainless steel                  | 4 047075 113101 |
| 4.676.40   | intertroch.plate for adolescents 90° 4.5 mm; 4 holes; arc depth 8 mm; blade 50 mm; stainless steel                   | 4 047075 143115 |
| <b>Plate intertroch. 90° 4.5 DC – plate thickness 5.4 mm</b>         |  |                 |
| 4.636.01   | intertroch.plate for adults with DC holes 90° 4.5 mm; 4 holes; arc depth 10 mm; blade 40 mm; stainless steel         | 4 047075 055500 |
| 4.636.02   | intertroch.plate for adults with DC holes 90° 4.5 mm; 4 holes; arc depth 10 mm; blade 50 mm; stainless steel         | 4 047075 055456 |
| 4.636.03   | intertroch.plate for adults with DC holes 90° 4.5 mm; 4 holes; arc depth 10 mm; blade 60 mm; stainless steel         | 4 047075 055524 |
| 4.636.04   | intertroch.plate for adults with DC holes 90° 4.5 mm; 4 holes; arc depth 15 mm; blade 40 mm; stainless steel         | 4 047075 055531 |
| 4.636.05   | intertroch.plate for adults with DC holes 90° 4.5 mm; 4 holes; arc depth 15 mm; blade 50 mm; stainless steel         | 4 047075 055548 |
| 4.636.06   | intertroch.plate for adults with DC holes 90° 4.5 mm; 4 holes; arc depth 15 mm; blade 60 mm; stainless steel         | 4 047075 110384 |
| 4.636.11   | intertroch.plate for adults with DC holes 90° 4.5 mm; 4 holes; arc depth 20 mm; blade 60 mm; stainless steel         | 4 047075 055395 |
| 4.636.29   | intertroch.plate for adults with DC holes 90° 4.5 mm; 4 holes; arc depth 10 mm; blade 45 mm; stainless steel         | 4 047075 055180 |
| <b>Plate intertroch. 90° 4.5/6.5 DC – plate thickness 5.4 mm</b>     |  |                 |
| 4.636.07SA   | intertroch.plate for adults with DC holes 90° 4.5 and 6.5 mm; 4 holes; arc depth 15 mm; blade 70 mm; stainless steel | 4 047075 113545 |

# Angle plates

stainless steel

| Item no.  | Item description  | EAN  |
|---|---|--|
| <b>Plate intertroch. 110° 4.5 DC – plate thickness 5.4 mm</b> |   |  |
| 4.656.00  | intertroch.plate for adults with DC holes 110° 4.5 mm; 4 holes; arc depth 10 mm; blade 60 mm; stainless steel | <br>4 047075 054718 |
| 4.656.01  | intertroch.plate for adults with DC holes 110° 4.5 mm; 4 holes; arc depth 10 mm; blade 65 mm; stainless steel | <br>4 047075 054541 |
| 4.656.015   | intertroch.plate for adults with DC holes 110° 4.5 mm; 4 holes; arc depth 10 mm; blade 50 mm; stainless steel | <br>4 047075 054558 |
| 4.656.02  | intertroch.plate for adults with DC holes 110° 4.5 mm; 4 holes; arc depth 10 mm; blade 75 mm; stainless steel | <br>4 047075 054660 |
| 4.656.022   | intertroch.plate for adults with DC holes 110° 4.5 mm; 4 holes; arc depth 10 mm; blade 80 mm; stainless steel | <br>4 047075 054657 |
| 4.656.025   | intertroch.plate for adults with DC holes 110° 4.5 mm; 4 holes; arc depth 10 mm; blade 55 mm; stainless steel | <br>4 047075 054589 |
| 4.656.027   | intertroch.plate for adults with DC holes 110° 4.5 mm; 4 holes; arc depth 10 mm; blade 70 mm; stainless steel | <br>4 047075 061884 |
| 4.656.03  | intertroch.plate for adults with DC holes 110° 4.5 mm; 4 holes; arc depth 10 mm; blade 85 mm; stainless steel | <br>4 047075 054671 |