


Set angled plates instruments

for angled plates; Set No. 19.230.00




Item no. 10.230.01



Item no. 10.230.05



Item no. 10.230.03



Item no. 19.230.00


Item no.	Item description	Quant.
Instruments		
10.230.01	triangular positioning plate 80°/70°/30°	1
10.230.02	triangular positioning plate 90°/50°/40°	1
10.230.03	triangular positioning plate 100°/60°/20°	1
10.230.04	triangular positioning plate for children 80°/70°/30°	1
10.230.05	triangular positioning plate for children 90°/50°/40°	1
10.230.06	triangular positioning plate for children 110°/60°/10°	1
10.230.07	positioning plate for variable osteotomy 110°/90°/90°/70°	1
10.230.08	positioning plate for variable osteotomy for children 110°/90°/90°/70°	1
10.230.09	condylar plate targeting device	1
10.230.10	impactor and extractor	1
10.230.11	impactor and extractor for angled plates for children	1
10.230.12	slit hammer	1
10.230.13	correcting bolt	1
10.230.14	plate fitting instrument for adults	1
10.230.15	plate fitting instrument for adolescents	1
10.230.16	plate fitting instrument for children	1
10.230.17	guide plate with adjustable angles	1
10.230.18	pin-countersink for three-jaw chuck	1
10.230.19	triple drill guide with fixed angle 130°	1
10.500.11	drill guide Ø3.2	1

Item no.	Item description	Quant.
Instruments		
10.500.20	drill guide for neutral and compression position, Ø3.2, inset 45	1
10.400.50	hex wrench, wrench width 5	1
2.904.08	drill bit for quick coupling, Ø3.2x145,- spiral length 60, double spiral	1
2.904.12	drill bit for quick coupling, Ø4.5x145,- spiral length 60, double spiral	1
2.925.45	tap for quick coupling for cortical screws, Ø4.5x130, thread 70	1
2.9406.35	hex screwdriver with handle for screws Ø4.5-Ø7.0, AF3.5, length 215	1
2.953.45	gauge with clasp for screws Ø4.5 and Ø6.5, measuring range 115	1
2.904.09	drill bit for quick coupling, Ø3.2x195,- spiral length 60, double spiral	1
2.904.13	drill bit for quick coupling, Ø4.5x195,- spiral length 85, double spiral	1
2.927.45	tap for quick coupling for cortical screws, Ø4.5x195, thread 70	1

Item no.	Item description	Quant.
Container		
19.231.00	perforated autoclavable container with inset for instruments for angled plates	1

Set angled plates instruments

for angled plates; Set No. 19.230.00



Item no. 19.230.00