

# Set of 2.4 vas distal radius plates

Variable angle-stable; for Ø 2.4 mm screws with TOBI® drive; titanium; set no. 19.117.20



Item no. 5.08173.04



Item no. 5.08193.04



Item no. 5.08273.04



Item no. 5.08293.04



Set no. 19.117.20 **TOBI**

Item no.	Item description	Quant.
<b>Implants - plates</b>		
<b>Variable angle-stable distal radius plates 2.4; titanium</b>		
5.08173.03	Distal radius plate 2.4; 7/3 holes; 23x55 mm; narrow; vas; left	1
5.08173.04	Distal radius plate 2.4; 7/4 holes; 23x63 mm; narrow; vas; left	1
5.08273.03	Distal radius plate 2.4; 7/3 holes; 23x55 mm; narrow; vas; right	1
5.08273.04	Distal radius plate 2.4; 7/4 holes; 23x63 mm; narrow; vas; right	1
5.08193.03	Distal radius plate 2.4; 9/3 holes; 28x55 mm; broad; vas; left	1
5.08193.04	Distal radius plate 2.4; 9/4 holes; 28x63 mm; broad; vas; left	1
5.08293.03	Distal radius plate 2.4; 9/3 holes; 28x55 mm; broad; vas; right	1
5.08293.04	Distal radius plate 2.4; 9/4 holes; 28x63 mm; broad; vas; right	1
<b>Instruments</b>		
2.904.33	Drill bit for QC; Ø1.8x105; Spiral length 30; double spiral	1
2.953.24	Gauge with clasp for Ø 2.4 screws; measuring range 60	1
2.9778.24	Angle-stable drill guide for drill 1.8 scaled; TOBI® drive B8; L 44; stainless steel	2
2.9286.08	Screwdriver for TOBI® drive B8 PPSU handle; L 187	1
10.500.24	Combi drill guide with handle for screwing 2.4, L 45	1
<b>Wires</b>		
6.031.15	Kirschner wire with trocar point and round end; Ø 1.5 mm x 150 mm; stainless steel	5
<b>Perforated container</b>		
19.118.20	Perforated container with insert for instruments and implants for 2.4 radius set	1

Item no.	Item description	Quant.
<b>Implants - screws</b>		
<b>Cortical screws Ø 2.4 mm with TOBI® drive; fully threaded; self tapping; titanium</b>		
3.2224.12	Length 12 mm	5
3.2224.14	Length 14 mm	5
3.2224.16	Length 16 mm	5
3.2224.18	Length 18 mm	5
3.2224.20	Length 20 mm	5

Item no.	Item description	Quant.
<b>Angle stable cortical screws Ø 2.4 mm with TOBI® drive; fully threaded; self tapping; titanium</b>		
3.2124.12	Length 12 mm	5
3.2124.14	Length 14 mm	5
3.2124.16	Length 16 mm	5
3.2124.18	Length 18 mm	5
3.2124.20	Length 20 mm	5
3.2124.22	Length 22 mm	5
3.2124.24	Length 24 mm	5
3.2124.26	Length 26 mm	5
3.2124.28	Length 28 mm	5

Item no.	Item description	Quant.
<b>Ø 1.8 mm angle stable pins with TOBI® drive; titanium</b>		
3.2324.16	Length 16 mm	3
3.2324.18	Length 18 mm	3
3.2324.20	Length 20 mm	3
3.2324.22	Length 22 mm	3
3.2324.24	Length 24 mm	3
3.2324.26	Length 26 mm	3
3.2324.28	Length 28 mm	3

# Set of 2.4 vas distal radius plates

Variable angle-stable; for Ø 2.4 mm screws with TOBI® drive; titanium; set no. 19.117.20



Set no. 19.117.20 **TOBI**

**Attention:** For the screws with TOBI® drive B8, only the corresponding screwdriver with the appropriate blade for TOBI® drive B8 (Item no. 2.9286.08) may be used.