


Set Distal fibula plate 3.5 vas IH T / IMPLANTS sterile

variable angle-stable; for screws Ø 3.5 mm; titanium

not sterile packed Set no. 19.474.00 / sterile packed Set no. 19.474.00PS+



Item no. 5.07803.09
Item no. S5.07803.09



All implants mentioned are available individually, not sterile and sterile packed.
Implants in the **Set not sterile packed** Set no. 19.474.00
Implants in **Set sterile packed** Set no. 19.474.00PS+

 **Internal hexagon**
Set no. 19.474.00
not sterile packed

Item no.	Item no. sterile	Item description	Quant.
Implantats – plates			
Fibula plate with conical thread 3.5, variable angle-stable, titanium			
5.07803.07	S5.07803.07	fibula plate with conical thread 3.5, 7 holes, 96 x 10, angle stable variable, ti	1
5.07803.09	S5.07803.09	fibula plate with conical thread 3.5, 9 holes, 122 x 10, angle stable variable, ti	1
5.07803.11	S5.07803.11	fibula plate with conical thread 3.5, 11 holes, 148 x 10, angle stable variable, ti	1
Implants – screws			
Cortical screw with conical head thread, Ø 3.5 mm, fully threaded, self tapping, titanium			 
3.133.10	S3.133.10	length 10	4
3.133.12	S3.133.12	length 12	4
3.133.14	S3.133.14	length 14	4
3.133.16	S3.133.16	length 16	4
3.133.18	S3.133.18	length 18	4
3.133.20	S3.133.20	length 20	4
3.133.22	S3.133.22	length 22	4
3.133.24	S3.133.24	length 24	4
3.133.26	S3.133.26	length 26	4
Cortical screw, Ø 3.5 mm, fully threaded, self tapping, titanium			 
3.132.10	S3.132.10	length 10	2
3.132.12	S3.132.12	length 12	2
3.132.14	S3.132.14	length 14	2
3.132.16	S3.132.16	length 16	2
3.132.18	S3.132.18	length 18	2
3.132.20	S3.132.20	length 20	2
3.132.22	S3.132.22	length 22	2
3.132.24	S3.132.24	length 24	2

Set Distal fibula plate 3.5 vas IH T / IMPLANTS sterile

variable angle-stable; for screws Ø 3.5 mm; titanium

not sterile packed Set no. 19.474.00 / sterile packed Set no. 19.474.00PS+

Item no.	Item no. sterile	Item description	Quant.
Implants – screws			
3.132.26	S3.132.26	length 26	2
3.132.28	S3.132.28	length 28	2
3.132.30	S3.132.30	length 30	2
3.132.34	S3.132.34	length 34	2
3.132.36	S3.132.36	length 36	2
3.132.38	S3.132.38	length 38	2

Instruments			
2.904.05	S2.904.05	drill bit for quick coupling, Ø 2.5 x 105 mm, spiral length 50 mm, double spiral	1
2.904.10	S2.904.10	drill bit for quick coupling, Ø 3.5 x 110 mm, spiral length 40 mm, double spiral	1
2.977.08	-	drill guide with handle for variable angle stable screwing, small fragment, length 45 mm	1
2.977.01	-	drill guide for angle stable screwing, length 45 mm, small fragment	2
10.500.27A	-	drill guide with handle for small fragment, handle length 55 mm	1
10.520.30	-	double drill guide Ø 2.5 mm/Ø 3.5 mm, spring bushing	1
2.9406.25	-	hex screwdriver with handle for screws Ø 2.7 mm - Ø 4.0 mm, AF 2.5, conical	1
2.953.60	-	gauge with clasp for screws Ø 3.5 mm, Ø 4.0 mm, measuring range 60 mm	1
2.911.06	-	bending lever for reco plates small fragment	2
2.954.01	-	screw forceps, self holding	1

Wires			
6.031.16	S6.031.16	Kirschner wire with trocar point and round end, Ø 1.6 x 150 mm, stainless steel	5

Container			
19.475.00/1	-	perforated autoclavable container with inset for instruments/implants angle stable plate small fragment	1

Set no. 19.474.00PS+ can be combined with:

19.474.00 Distal Fibula plate 3.5 vas IH T SET
19.035.00I SF 3.5 tools IH INSTRUMENTS



Königsee Implantate GmbH
OT Aschau | Am Sand 4
07426 Allendorf | Germany


Fon +49 36738 498-560
Fax +49 36738 498-569

Mail export@koenigsee-implantate.de

For further information please visit:
www.koenigsee-implantate.de

StZi | Edition 03 | 2024-07
All pictured products are CE compliant.

Set Distal fibula plate 3.5 vas IH T / IMPLANTS sterile
variable angle-stable; for screws Ø 3.5 mm; titanium
not sterile packed Set no. 19.474.00 / sterile packed Set no. 19.474.00PS+



 Internal hexagon
Set no. 19.474.00
not sterile packed