


# Set VaWiKo® distal fibula plate 3.5 IH T / IMPLANTS sterile

variable angle-stable; for screws Ø 3.5 mm; head 5.0 mm; titanium


not sterile packed Set no. 19.474.00/2 / sterile packed Set no. 19.474.00/2PS+







Item no. 5.07813.09  
Item no. S5.07813.09



All implants mentioned are available individually, not sterile and sterile packed.  
Implants in the **Set not sterile packed** Set no. 19.474.00/2  
Implants in **Set sterile packed** Set no. 19.474.00/2PS+

 **Internal hexagon**  
Set no. 19.474.00/2  
not sterile packed

| Item no.  | Item no. sterile | Item description   | Quant.  |
|---|------------------|--|---|
| <b>Implantats – plates</b>  |                  |  |   |
| <b>VaWiKo® distal fibula plate 3.5 holeØ5; titanium</b>   |                  |  |   |
| 5.07813.06  | S5.07813.06      | VaWiKo® lat. dist. Fibula Pl. 3.5; 4/2 hole; length 83 mm  | 1   |
| 5.07813.07  | S5.07813.07      | VaWiKo® lat. dist. Fibula Pl. 3.5; 4/3 hole; length 96 mm  | 1   |
| 5.07813.09  | S5.07813.09      | VaWiKo® lat. dist. Fibula Pl. 3.5; 4/5 hole; length 122 mm | 1   |
| <b>Implants – screws</b>  |                  |  |   |
| <b>Cortical screw with conical head thread Ø 3.5 mm; fully threaded; self tapping; head 5.0; titanium</b> |                  |  |   |
| 3.136.10  | S3.136.10        | length 10  | 4   |
| 3.136.12  | S3.136.12        | length 12  | 4   |
| 3.136.14  | S3.136.14        | length 14  | 4   |
| 3.136.16  | S3.136.16        | length 16  | 4   |
| 3.136.18  | S3.136.18        | length 18  | 4   |
| 3.136.20  | S3.136.20        | length 20  | 4   |
| 3.136.22  | S3.136.22        | length 22  | 4   |
| 3.136.24  | S3.136.24        | length 24  | 4   |
| <b>Cortical screw Ø 3.5 mm; fully threaded; self tapping, titanium</b>                                    |                  |  |   |
| 3.138.10  | S3.138.10        | length 10  | 2   |
| 3.138.12  | S3.138.12        | length 12  | 2   |
| 3.138.14  | S3.138.14        | length 14  | 2   |
| 3.138.16  | S3.138.16        | length 16  | 2   |
| 3.138.18  | S3.138.18        | length 18  | 2   |
| 3.138.20  | S3.138.20        | length 20  | 2   |
| 3.138.22  | S3.138.22        | length 22  | 2   |
| 3.138.24  | S3.138.24        | length 24  | 2   |

# Set VaWiKo® distal fibula plate 3.5 IH T / IMPLANTS sterile

variable angle-stable; for screws Ø 3.5 mm; head 5.0 mm; titanium

not sterile packed Set no. 19.474.00/2 / sterile packed Set no. 19.474.00/2PS+

| Item no.          | Item no. sterile | Item description | Quant. |
|-------------------|------------------|------------------|--------|
| Implants – screws |                  |                  |        |
| 3.138.36          | S3.138.36        | length 36        | 2      |
| 3.138.38          | S3.138.38        | length 38        | 2      |
| 3.138.40          | S3.138.40        | length 40        | 2      |
| 3.138.42          | S3.138.42        | length 42        | 2      |
| 3.138.44          | S3.138.44        | length 44        | 2      |
| 3.138.46          | S3.138.46        | length 46        | 2      |
| 3.138.48          | S3.138.48        | length 48        | 2      |
| 3.138.50          | S3.138.50        | length 50        | 2      |

|             |           |   |   |
|-------------|-----------|---|---|
| Instruments |           |   |   |
| 2.904.06    | S2.904.06 | drill bit for quick coupling, Ø 2.5 x 125 mm, spiral length 60 mm, double spiral  | 1 |
| 2.904.10    | S2.904.10 | drill bit for quick coupling, Ø 3.5 x 110 mm, spiral length 40 mm, double spiral  | 1 |
| 2.977.12    | -         | drill guide with handle for variable angle stable screwing 3.5 with head Ø 5.0 mm | 1 |
| 2.977.13    | -         | drill guide for angle stable screws 3.5 with head Ø 5.0 mm, length 44 mm          | 2 |
| 10.520.30   | -         | double drill guide Ø 2.5 mm/ Ø 3.5 mm, spring bushing                             | 1 |
| 2.9406.25   | -         | hex screwdriver with handle for screws Ø 2.7 mm - Ø 4.0 mm, AF 2.5, conical       | 1 |
| 2.953.60    | -         | gauge with clasp for screws Ø 3.5 mm, Ø 4.0 mm, measuring range 60 mm             | 1 |
| 2.924.35    | -         | tap for quick coupling for cortical screws, Ø 3.5 x 110 mm                        | 1 |
| 2.954.01    | -         | screw forceps, self holding   | 1 |

|             |           |  |   |
|-------------|-----------|--|---|
| Accessories |           |  |   |
| 3.200.02    | S3.200.02 | washer, outside Ø 7.0 mm / inside Ø 3.8 mm, titanium | 3 |

|          |           |   |   |
|----------|-----------|---|---|
| Wires    |           |   |   |
| 6.031.16 | S6.031.16 | Kirschner wire with trocar point and round end, Ø 1.6 x 150 mm, stainless steel | 5 |

|             |   |   |   |
|-------------|---|---|---|
| Container   |   |   |   |
| 19.475.00/2 | - | perforated autoclav. container with inset for instruments/ implants vas, plate and screws up to 50 mm in length | 1 |

#### Optionally available:

5.07813.11 S5.07813.11 VaWiKo® lat. dist. Fibula Pl. 3.5; 4/7 hole; length 148 mm

#### Set no. 19.474.00/2PS+ can be combined with:

19.474.00/2 VaWiKo®Distal Fibula plate 3.5 IH T SET  
19.035.001 SF 3.5 tools IH INSTRUMENTS



Königsee Implantate GmbH  
OT Aschau | Am Sand 4  
07426 Allendorf | Germany

Fon +49 36738 498-560  
Fax +49 36738 498-569

Mail [export@koenigsee-implantate.de](mailto:export@koenigsee-implantate.de)


For further information please visit:  
[www.koenigsee-implantate.de](http://www.koenigsee-implantate.de)


StZi | Edition 04 | 2024-07  
All pictured products are CE compliant.

# Set VaWiKo® distal fibula plate 3.5 IH T / IMPLANTS sterile

variable angle-stable; for screws  $\varnothing$  3.5 mm; head 5.0 mm; titanium

not sterile packed Set no. 19.474.00/2 / sterile packed Set no. 19.474.00/2PS+



 Internal hexagon  
Set no. 19.474.00/2  
not sterile packed