

Cancellous screws

angle stable; titanium

not sterile packed / sterile packed



Item no. 3.147.40
Item no. S3.147.40

Item no.	EAN	Description	Item no. sterile	EAN sterile
Cancellous screw 4.0 as st				
3.147.10		Cancellous screw; diameter 4.0 mm x 10 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.10	
3.147.12		Cancellous screw; diameter 4.0 mm x 12 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.12	
3.147.14		Cancellous screw; diameter 4.0 mm x 14 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.14	
3.147.16		Cancellous screw; diameter 4.0 mm x 16 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.16	
3.147.18		Cancellous screw; diameter 4.0 mm x 18 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.18	
3.147.20		Cancellous screw; diameter 4.0 mm x 20 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.20	
3.147.22		Cancellous screw; diameter 4.0 mm x 22 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.22	
3.147.24		Cancellous screw; diameter 4.0 mm x 24 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.24	
3.147.26		Cancellous screw; diameter 4.0 mm x 26 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.26	
3.147.28		Cancellous screw; diameter 4.0 mm x 28 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.28	
3.147.30		Cancellous screw; diameter 4.0 mm x 30 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.30	
3.147.32		Cancellous screw; diameter 4.0 mm x 32 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.32	
3.147.34		Cancellous screw; diameter 4.0 mm x 34 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.34	
3.147.35		Cancellous screw; diameter 4.0 mm x 35 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.35	
3.147.36		Cancellous screw; diameter 4.0 mm x 36 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.36	
3.147.38		Cancellous screw; diameter 4.0 mm x 38 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.38	
3.147.40		Cancellous screw; diameter 4.0 mm x 40 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.40	

Cancellous screws

angle stable; titanium

not sterile packed / sterile packed



Item no. 3.147.40
Item no. S3.147.40

Item no.	EAN	Description	Item no. sterile	EAN sterile
Cancellous screw 4.0 as st				
3.147.42		Cancellous screw; diameter 4.0 mm x 42 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.42	
3.147.44		Cancellous screw; diameter 4.0 mm x 44 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.44	
3.147.45		Cancellous screw; diameter 4.0 mm x 45 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.45	
3.147.46		Cancellous screw; diameter 4.0 mm x 46 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.46	
3.147.48		Cancellous screw; diameter 4.0 mm x 48 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.48	
3.147.50		Cancellous screw; diameter 4.0 mm x 50 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.50	
3.147.52		Cancellous screw; diameter 4.0 mm x 52 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.52	
3.147.54		Cancellous screw; diameter 4.0 mm x 54 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.54	
3.147.55		Cancellous screw; diameter 4.0 mm x 55 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.55	
3.147.56		Cancellous screw; diameter 4.0 mm x 56 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.56	
3.147.58		Cancellous screw; diameter 4.0 mm x 58 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.58	
3.147.60		Cancellous screw; diameter 4.0 mm x 60 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.60	
3.147.65		Cancellous screw; diameter 4.0 mm x 65 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.65	
3.147.70		Cancellous screw; diameter 4.0 mm x 70 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.70	
3.147.75		Cancellous screw; diameter 4.0 mm x 75 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.75	
3.147.80		Cancellous screw; diameter 4.0 mm x 80 mm; angle stable; fully threaded; self tapping; internal hexagon 2.5 mm; titanium	S3.147.80	

Cancellous screws

angle stable; titanium

not sterile packed / sterile packed



Item no. 3.0340.35
Item no. S3.0340.35

Item no.	EAN	Description	Item no. sterile	EAN sterile
Cancellous screw ² 4.0 as SL15 st				
3.0340.10		Cancellous screw ² ; diameter 4.0 mm x 10 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.10	
3.0340.12		Cancellous screw ² ; diameter 4.0 mm x 12 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.12	
3.0340.14		Cancellous screw ² ; diameter 4.0 mm x 14 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.14	
3.0340.16		Cancellous screw ² ; diameter 4.0 mm x 16 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.16	
3.0340.18		Cancellous screw ² ; diameter 4.0 mm x 18 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.18	
3.0340.20		Cancellous screw ² ; diameter 4.0 mm x 20 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.20	
3.0340.22		Cancellous screw ² ; diameter 4.0 mm x 22 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.22	
3.0340.24		Cancellous screw ² ; diameter 4.0 mm x 24 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.24	
3.0340.26		Cancellous screw ² ; diameter 4.0 mm x 26 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.26	
3.0340.28		Cancellous screw ² ; diameter 4.0 mm x 28 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.28	
3.0340.30		Cancellous screw ² ; diameter 4.0 mm x 30 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.30	
3.0340.32		Cancellous screw ² ; diameter 4.0 mm x 32 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.32	
3.0340.34		Cancellous screw ² ; diameter 4.0 mm x 34 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.34	

Cancellous screws

angle stable; titanium

not sterile packed / sterile packed



Item no. 3.0340.35
Item no. S3.0340.35

Item no.	EAN	Description	Item no. sterile	EAN sterile
Cancellous screw ² 4.0 as SL15 st				
3.0340.35		Cancellous screw ² ; diameter 4.0 mm x 35 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.35	
3.0340.36		Cancellous screw ² ; diameter 4.0 mm x 36 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.36	
3.0340.38		Cancellous screw ² ; diameter 4.0 mm x 38 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.38	
3.0340.40		Cancellous screw ² ; diameter 4.0 mm x 40 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.40	
3.0340.42		Cancellous screw ² ; diameter 4.0 mm x 42 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.42	
3.0340.44		Cancellous screw ² ; diameter 4.0 mm x 44 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.44	
3.0340.45		Cancellous screw ² ; diameter 4.0 mm x 45 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.45	
3.0340.46		Cancellous screw ² ; diameter 4.0 mm x 46 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.46	
3.0340.48		Cancellous screw ² ; diameter 4.0 mm x 48 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.48	
3.0340.50		Cancellous screw ² ; diameter 4.0 mm x 50 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.50	
3.0340.52		Cancellous screw ² ; diameter 4.0 mm x 52 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.52	
3.0340.54		Cancellous screw ² ; diameter 4.0 mm x 54 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.54	
3.0340.55		Cancellous screw ² ; diameter 4.0 mm x 55 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.55	

Cancellous screws

angle stable; titanium

not sterile packed / sterile packed



Item no. 3.0340.35
Item no. S3.0340.35



Item no. 3.163.45
Item no. S3.163.45

Item no.	EAN	Description	Item no. sterile	EAN sterile
Cancellous screw² 4.0 as SL15 st				
3.0340.56		Cancellous screw ² ; diameter 4.0 mm x 56 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.56	
3.0340.58		Cancellous screw ² ; diameter 4.0 mm x 58 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.58	
3.0340.60		Cancellous screw ² ; diameter 4.0 mm x 60 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.60	
3.0340.65		Cancellous screw ² ; diameter 4.0 mm x 65 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.65	
3.0340.70		Cancellous screw ² ; diameter 4.0 mm x 70 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.70	
3.0340.75		Cancellous screw ² ; diameter 4.0 mm x 75 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.75	
3.0340.80		Cancellous screw ² ; diameter 4.0 mm x 80 mm; angle stable; fully threaded; self tapping; double threaded head; six lobe 15; titanium	S3.0340.80	
Cancellous screw 6.0 as st				
3.163.24		Cancellous screw; diameter 6.0 mm x 24 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.24	
3.163.26		Cancellous screw; diameter 6.0 mm x 26 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.26	
3.163.28		Cancellous screw; diameter 6.0 mm x 28 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.28	
3.163.30		Cancellous screw; diameter 6.0 mm x 30 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.30	
3.163.32		Cancellous screw; diameter 6.0 mm x 32 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.32	
3.163.34		Cancellous screw; diameter 6.0 mm x 34 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.34	
3.163.35		Cancellous screw; diameter 6.0 mm x 35 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.35	

Cancellous screws

angle stable; titanium

not sterile packed / sterile packed



Item no. 3.163.45
Item no. S3.163.45

Item no.	EAN	Description	Item no. sterile	EAN sterile
Cancellous screw 6.0 as st				
3.163.36		Cancellous screw; diameter 6.0 mm x 36 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.36	
3.163.38		Cancellous screw; diameter 6.0 mm x 38 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.38	
3.163.40		Cancellous screw; diameter 6.0 mm x 40 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.40	
3.163.42		Cancellous screw; diameter 6.0 mm x 42 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.42	
3.163.44		Cancellous screw; diameter 6.0 mm x 44 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.44	
3.163.45		Cancellous screw; diameter 6.0 mm x 45 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.45	
3.163.46		Cancellous screw; diameter 6.0 mm x 46 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.46	
3.163.48		Cancellous screw; diameter 6.0 mm x 48 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.48	
3.163.50		Cancellous screw; diameter 6.0 mm x 50 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.50	
3.163.52		Cancellous screw; diameter 6.0 mm x 52 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.52	
3.163.54		Cancellous screw; diameter 6.0 mm x 54 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.54	
3.163.55		Cancellous screw; diameter 6.0 mm x 55 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.55	
3.163.56		Cancellous screw; diameter 6.0 mm x 56 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.56	
3.163.58		Cancellous screw; diameter 6.0 mm x 58 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.58	
3.163.60		Cancellous screw; diameter 6.0 mm x 60 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.60	
3.163.62		Cancellous screw; diameter 6.0 mm x 62 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.62	
3.163.64		Cancellous screw; diameter 6.0 mm x 64 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.64	
3.163.65		Cancellous screw; diameter 6.0 mm x 65 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.65	

Cancellous screws

angle stable; titanium

not sterile packed / sterile packed



Item no. 3.163.45
Item no. S3.163.45

Item no.	EAN	Description	Item no. sterile	EAN sterile
Cancellous screw 6.0 as st				
3.163.66		Cancellous screw; diameter 6.0 mm x 66 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.66	
3.163.68		Cancellous screw; diameter 6.0 mm x 68 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.68	
3.163.70		Cancellous screw; diameter 6.0 mm x 70 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.70	
3.163.72		Cancellous screw; diameter 6.0 mm x 72 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.72	
3.163.74		Cancellous screw; diameter 6.0 mm x 74 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.74	
3.163.75		Cancellous screw; diameter 6.0 mm x 75 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.75	
3.163.76		Cancellous screw; diameter 6.0 mm x 76 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.76	
3.163.78		Cancellous screw; diameter 6.0 mm x 78 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.78	
3.163.80		Cancellous screw; diameter 6.0 mm x 80 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.80	
3.163.82		Cancellous screw; diameter 6.0 mm x 82 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.82	
3.163.84		Cancellous screw; diameter 6.0 mm x 84 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.84	
3.163.85		Cancellous screw; diameter 6.0 mm x 85 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.85	
3.163.86		Cancellous screw; diameter 6.0 mm x 86 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.86	
3.163.90		Cancellous screw; diameter 6.0 mm x 90 mm; angle stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium	S3.163.90	