


Set Prox.Humerus EIGHT² vas Corti SL T / IMPLANTS sterile

double threaded; variable angle-stable; for screws Ø 4.0/3.5 mm; titanium

not sterile packed Set no. 19.3683.10 / sterile packed Set no. 19.3683.10PS+



Item no. 5.07933.86
Item no. S5.07933.86



All implants mentioned are available individually, not sterile and sterile packed.
Implants in the **Set not sterile packed** Set no. 19.3683.10
Implants in **Set sterile packed** Set no. 19.3683.10PS+

★ Six lobe
Set no. 19.3683.10
not sterile packed

Item no.	Item no. sterile	Item description	Quant.
Implantats – plates			
Proximal humerus plate ² 4.0/3.5			
5.07933.83	S5.07933.83	prox. humerus plate 4.0/3.5; 8/3H; L 77 mm; vas; 2t; T	1
5.07933.84	S5.07933.84	prox. humerus plate 4.0/3.5; 8/4H; L 93 mm; vas; 2t; T	1
5.07933.85	S5.07933.85	prox. humerus plate 4.0/3.5; 8/5H; L 109 mm; vas; 2t; T	1
5.07933.86	S5.07933.86	prox. humerus plate 4.0/3.5; 8/6H; L 125 mm; vas; 2t; T	1
Implants – screws			
Cortical screw 3.5 SL15 st  ★			
3.1321.22	S3.1321.22	Cortical screw; Ø 3.5 x 22; FT; SL15; st; T	5
3.1321.24	S3.1321.24	Cortical screw; Ø 3.5 x 24; FT; SL15; st; T	5
3.1321.26	S3.1321.26	Cortical screw; Ø 3.5 x 26; FT; SL15; st; T	5
3.1321.28	S3.1321.28	Cortical screw; Ø 3.5 x 26; FT; SL15; st; T	5
3.1321.30	S3.1321.30	Cortical screw; Ø 3.5 x 30; FT; SL15; st; T	5
3.1321.32	S3.1321.32	Cortical screw; Ø 3.5 x 32; FT; SL15; st; T	5
3.1321.34	S3.1321.34	Cortical screw; Ø 3.5 x 34; FT; SL15; st; T	5
3.1321.36	S3.1321.36	Cortical screw; Ø 3.5 x 36; FT; SL15; st; T	5
3.1321.38	S3.1321.38	Cortical screw; Ø 3.5 x 38; FT; SL15; st; T	5
3.1321.40	S3.1321.40	Cortical screw; Ø 3.5 x 40; FT; SL15; st; T	5
3.1321.42	S3.1321.42	Cortical screw; Ø 3.5 x 42; FT; SL15; st; T	5
Cortical screw ² 3.5 as SL15 st  ★			
3.0135.22	S3.0135.22	Cortical screw ² ; as; Ø 3.5 x 22; FT; SL15; st; T	5
3.0135.24	S3.0135.24	Cortical screw ² ; as; Ø 3.5 x 24; FT; SL15; st; T	5
3.0135.26	S3.0135.26	Cortical screw ² ; as; Ø 3.5 x 26; FT; SL15; st; T	5
3.0135.28	S3.0135.28	Cortical screw ² ; as; Ø 3.5 x 28; FT; SL15; st; T	5
3.0135.30	S3.0135.30	Cortical screw ² ; as; Ø 3.5 x 30; FT; SL15; st; T	5

Set Prox.Humerus EIGHT² vas Corti SL T / IMPLANTS sterile

double threaded; variable angle-stable; for screws Ø 4.0/3.5 mm; titanium

not sterile packed Set no. 19.3683.10 / sterile packed Set no. 19.3683.10PS+


Item no.	Item no. sterile	Item description	Quant.
Implants – screws			
3.0135.32	S3.0135.32	Cortical screw ² ; as; Ø 3.5 x 32; FT; SL15; st; T	5
3.0135.34	S3.0135.34	Cortical screw ² ; as; Ø 3.5 x 34; FT; SL15; st; T	5
3.0135.36	S3.0135.36	Cortical screw ² ; as; Ø 3.5 x 36; FT; SL15; st; T	5
3.0135.38	S3.0135.38	Cortical screw ² ; as; Ø 3.5 x 38; FT; SL15; st; T	5
3.0135.40	S3.0135.40	Cortical screw ² ; as; Ø 3.5 x 40; FT; SL15; st; T	5
3.0135.42	S3.0135.42	Cortical screw ² ; as; Ø 3.5 x 42; FT; SL15; st; T	5
Cortical screw² 4.0 as SL15 st 			
3.0140.34	S3.0140.34	cortical screw; Ø 4.0 x 34 mm; angle-stable; fully threaded; s.t.	5
3.0140.36	S3.0140.36	cortical screw; Ø 4.0 x 36 mm; angle-stable; fully threaded; s.t.	5
3.0140.38	S3.0140.38	cortical screw; Ø 4.0 x 38 mm; angle-stable; fully threaded; s.t.	5
3.0140.40	S3.0140.40	cortical screw; Ø 4.0 x 40 mm; angle-stable; fully threaded; s.t.	5
3.0140.42	S3.0140.42	cortical screw; Ø 4.0 x 42 mm; angle-stable; fully threaded; s.t.	5
3.0140.44	S3.0140.44	cortical screw; Ø 4.0 x 44 mm; angle-stable; fully threaded; s.t.	5
3.0140.46	S3.0140.46	cortical screw; Ø 4.0 x 46 mm; angle-stable; fully threaded; s.t.	5
3.0140.48	S3.0140.48	cortical screw; Ø 4.0 x 48 mm; angle-stable; fully threaded; s.t.	5
3.0140.50	S3.0140.50	cortical screw; Ø 4.0 x 50 mm; angle-stable; fully threaded; s.t.	5
3.0140.52	S3.0140.52	cortical screw; Ø 4.0 x 52 mm; angle-stable; fully threaded; s.t.	5
3.0140.54	S3.0140.54	cortical screw; Ø 4.0 x 54 mm; angle-stable; fully threaded; s.t.	5
3.0140.56	S3.0140.56	cortical screw; Ø 4.0 x 56 mm; angle-stable; fully threaded; s.t.	5
Instruments			
10.173.05	—	screw driver for guide block	1
10.179.11	—	Aiming block humerus plate EIGHT 5.07933.XX; Ln 54; PEEK	1
10.173.01	—	sleeve for guide block for proximal small fragment humerus plates	2
2.904.22	—	Drill bit qc 2.5; Spiral 60; 2-flute; L 155; green	1
2.904.25	—	Drill bit qc 3.0; Spiral 60; 2-flute; L 155	1
10.500.27A	—	drill guide with handle for small fragment, handle length 55	1
2.9385.15	—	screwdriver with handle SL15; Ln 215 mm	1
2.953.70	—	gauge w. clasp for screws Ø 3.5, Ø 4.0, measuring range 70	1
2.954.01	—	screw forceps, self holdingte	1
2.977.23	—	Drill guide as for drill 2.5 double threaded, Ln 68, StainlessSteel	1
2.977.09	—	drill guide w.handle for variable angle stable screwing, small fragment, length 68	1
2.977.22	—	drill guide angle stable screwing 4.0; L 68 mm; 2 t	1
2.977.21	—	drill guide with handle vas screwing 4.0; L 68mm	1
2.9776.30	—	sinkhole drill guide with handle; drill Ø 3.0/ Ø 2.5	1

Set Prox.Humerus EIGHT² vas Corti SL T / IMPLANTS sterile

double threaded; variable angle-stable; for screws Ø 4.0/3.5 mm; titanium

not sterile packed Set no. 19.3683.10 / sterile packed Set no. 19.3683.10PS+

Item no.	Item no. sterile	Item description	Quant.
Wires			
6.031.18	—	Kirschner wire with trocar point and round end, Ø 1.8 x 150, st.st.	5
Container			
19.3661.10	—	Perforated container humerus set double threaded	1



★ Six lobe
Set no. 19.3683.10
not sterile packed