

Set Screw container for DHS

stainless steel; set no. 19.168.10




Item no. 2.152.30



Item no. 2.170.50



Item no. 2.180.50




Set no. 19.168.10 *  Internal hexagon

Item no.	Item description	Quant.
Implants – screws		
Cortical screw Ø 4.5 mm; fully threaded; self tapping; stainless steel		
2.152.20	length 20 mm	3
2.152.22	length 22 mm	3
2.152.24	length 24 mm	3
2.152.26	length 26 mm	3
2.152.28	length 28 mm	3
2.152.30	length 30 mm	3
2.152.32	length 32 mm	3
2.152.34	length 34 mm	3
2.152.36	length 36 mm	3
2.152.38	length 38 mm	3
2.152.40	length 40 mm	3
2.152.42	length 42 mm	3
2.152.44	length 44 mm	3
2.152.46	length 46 mm	3
2.152.48	length 48 mm	3
2.152.50	length 50 mm	2
2.152.52	length 52 mm	2
2.152.54	length 54 mm	2
2.152.56	length 56 mm	2
2.152.58	length 58 mm	2
2.152.60	length 60 mm	2

Item no.	Item description	Quant.
Implants – screws		
Cancellous screw Ø 6.5 mm; thread 16; stainless steel		
2.170.45	length 45 mm	2
2.170.50	length 50 mm	2
2.170.55	length 55 mm	2
2.170.60	length 60 mm	2
2.170.65	length 65 mm	2
2.170.70	length 70 mm	2
2.170.75	length 75 mm	2
2.170.80	length 80 mm	2
2.170.85	length 85 mm	2
2.170.90	length 90 mm	2
2.170.95	length 95 mm	2
2.170.100	length 100 mm	2
2.170.105	length 105 mm	2
2.170.110	length 110 mm	2

Set Screw container for DHS

stainless steel; set no. 19.168.10



Set no. 19.168.10 *  Internal hexagon

Item no.	Item description	Quant.
Implants – screws		
Cancellous screw Ø 6.5 mm; thread 32; stainless steel		
2.180.45	length 45 mm	2
2.180.50	length 50 mm	2
2.180.55	length 55 mm	2
2.180.60	length 60 mm	2
2.180.65	length 65 mm	2
2.180.70	length 70 mm	2
2.180.75	length 75 mm	2
2.180.80	length 80 mm	2
2.180.85	length 85 mm	2
2.180.90	length 90 mm	2
2.180.95	length 95 mm	2
2.180.100	length 100 mm	2
2.180.105	length 105 mm	2
2.180.110	length 110 mm	2

Item no.	Item description	Quant.
Instruments		
2.954.01	screw forceps; self holding	1
20.502.00	holding clamp; big	1

Implants – washers		
2.200.03	washer; outside Ø 13 mm / inside Ø 6.7 mm; stainless steel	3

Container		
19.168.00	screw container for large fragment; empty	1

* Example illustration can show different number of screws.