

Set Screw Container for DHS

titanium; Set no. 19.168.30




Item no. 3.152.30



Item no. 3.170.70



Item no. 3.180.70




Set no. 19.168.30 *  internal hexagon

Item no.	Item description	Quant.
Implants – screws		
Cortical screw Ø 4.5 mm; fully threaded; self tapping; titanium		
3.152.20	length 20 mm	3
3.152.22	length 22 mm	3
3.152.24	length 24 mm	3
3.152.26	length 26 mm	3
3.152.28	length 28 mm	3
3.152.30	length 30 mm	3
3.152.32	length 32 mm	3
3.152.34	length 34 mm	3
3.152.36	length 36 mm	3
3.152.38	length 38 mm	3
3.152.40	length 40 mm	3
3.152.42	length 42 mm	3
3.152.44	length 44 mm	3
3.152.46	length 46 mm	3
3.152.48	length 48 mm	3
3.152.50	length 50 mm	2
3.152.52	length 52 mm	2
3.152.54	length 54 mm	2
3.152.56	length 56 mm	2
3.152.58	length 58 mm	2
3.152.60	length 60 mm	2

Item no.	Item description	Quant.
Implants – screws		
Cancellous screw Ø 6.5 mm; thread 16; titanium		
3.170.45	length 45 mm	2
3.170.50	length 50 mm	2
3.170.55	length 55 mm	2
3.170.60	length 60 mm	2
3.170.65	length 65 mm	2
3.170.70	length 70 mm	2
3.170.75	length 75 mm	2
3.170.80	length 80 mm	2
3.170.85	length 85 mm	2
3.170.90	length 90 mm	2
3.170.95	length 95 mm	2
3.170.100	length 100 mm	2
3.170.105	length 105 mm	2
3.170.110	length 110 mm	2

Set Screw Container for DHS

titanium; Set no. 19.168.30



Set no. 19.168.30 *  internal hexagon

Item no.	Item description	Quant.
Implants – screws		
Cancellous screw Ø 6.5 mm; thread 32; titanium		
3.180.45	length 45 mm	2
3.180.50	length 50 mm	2
3.180.55	length 55 mm	2
3.180.60	length 60 mm	2
3.180.65	length 65 mm	2
3.180.70	length 70 mm	2
3.180.75	length 75 mm	2
3.180.80	length 80 mm	2
3.180.85	length 85 mm	2
3.180.90	length 90 mm	2
3.180.95	length 95 mm	2
3.180.100	length 100 mm	2
3.180.105	length 105 mm	2
3.180.110	length 110 mm	2

Item no.	Item description	Quant.
Instruments		
2.954.01	screw forceps, self holding	1
20.502.00	holding clamp, big	1

Implants – washers		
3.200.03	washer, outside Ø 13 mm / inside Ø 6.7 mm, titanium	3

Container		
19.168.00	screw container for large fragment; empty	1

* Example illustration can show different number of screws.