

Set Instrumental Tornillos canulados 3.5 HI

para tornillos canulados Ø 3.5; N.º set 19.035.201



Ilustración 19.035.201  Hexágono interior

N.º art.	Descripción	Cant.
Instrumental		
10.935.00	Sistema d.casquillo taladro percután.p.torn.perforados Ø3,5	1
2.905.27	Broca helicoidal ACRA2.7; espiral30; labio4; L160; can.1.35	1
2.905.35	Broca helicoidal ACRA3.5; espiral30; labio4; L160; can.1.35	1
2.907.35	Fresa p.espacio d.cab.; Ø5x160; perf.,lØ1,3; p.torn.perfor. Ø3,5	1
10.306.06	Empuñadura en T c.acoplamiento rápido; L 93,5; perforada,lØ2,5	1
2.9436.25	Dest. hexagonal c.mango p.Ø2,7 a Ø4,0; perf.; lØ1.3; AnLi2.5; L 215mm	1
2.954.27A	Casquillo de sujeción para tornillos Ø2,7 y Ø3,5	1
2.950.150	Varilla de medición p. torn. perforados Ø2,7-Ø4,0; L 150; camp. med. 80	1

N.º art.	Descripción	Cant.
Alambres		
2.978.12	Aguja guía;Ø1,2x150;rosca 5;A	5

Tambor perforado		
19.040.01	Tambor perf. para Instrumental tornillos can. 2.7-4.5	1

Set Instrumental Tornillos canulados 3.5 HI

para tornillos canulados Ø 3.5; N.º set 19.035.20I

Puede combinarse con:

19.184.00	Tornillos can. 3.5 ar ap HI T SET
19.184.00PS+	Tornillos can. 3.5 ar ap HI T TORNILLOS estéril
19.181.00	Tornillos can. 3.5 ar HI T SET
19.181.00PS+	Tornillos can. 3.5 ar HI T TORNILLOS estéril




Ilustración 19.035.20I  Hexágono interior