

# Set Femur distal groove and notch IH S / IMPLANTS sterile

angle-stable; for screws Ø 4.5/5.5 mm; stainless steel


not sterile packed Set no. 19.600.10 / sterile packed Set no. 19.600.10PS+



Item no. 4.9502.12  
Item no. S4.9502.12




Item no. 4.9512.12  
Item no. S4.9512.12




All implants mentioned are available individually, not sterile and sterile packed.

Implants in the set **not sterile packed** Set no. 19.600.10

Implants in the set **sterile packed** Set no. 19.600.10PS+

 Internal hexagon  
Set no. 19.600.10  
not sterile packed



Item no.	Item no. sterile	Item description	Quant.
<b>Implants – plates</b>			
<b>Distal femur plate 4.5-6.0 as notch</b>			
4.9502.12	S4.9502.12	Distal femur plate with grooves and notches; 12 holes; 246 mm x 18 mm; right; stainless steel	1
4.9502.16	S4.9502.16	Distal femur plate with grooves and notches; 16 holes; 309 mm x 18 mm; right; stainless steel	1
4.9512.12	S4.9512.12	Distal femur plate with grooves and notches; 12 holes; 246 mm x 18 mm; left; stainless steel	1
4.9512.16	S4.9512.16	Distal femur plate with grooves and notches; 16 holes; 309 mm x 18 mm; left; stainless steel	1

<b>Implants – screws</b>			
<b>Cortical screw 4.5 st</b> 			
2.152.14	S2.152.14	Ø 4.5 mm x 14 mm; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.152.18	S2.152.18	Ø 4.5 mm x 18 mm; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.152.22	S2.152.22	Ø 4.5 mm x 22 mm; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.152.26	S2.152.26	Ø 4.5 mm x 26 mm; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.152.30	S2.152.30	Ø 4.5 mm x 30 mm; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.152.34	S2.152.34	Ø 4.5 mm x 34 mm; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.152.38	S2.152.38	Ø 4.5 mm x 38 mm; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.152.42	S2.152.42	Ø 4.5 mm x 42 mm; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.152.46	S2.152.46	Ø 4.5 mm x 46 mm; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.152.50	S2.152.50	Ø 4.5 mm x 50 mm; fully threaded; self tapping; internal hexagon 3.5; stainless steel	3
2.152.54	S2.152.54	Ø 4.5 mm x 54 mm; fully threaded; self tapping; internal hexagon 3.5; stainless steel	3

# Set Femur distal groove and notch IH S / IMPLANTS sterile

angle-stable; for screws Ø 4.5/5.5 mm; stainless steel



not sterile packed Set no. 19.600.10 / sterile packed Set no. 19.600.10PS+

Item no.	Item no. sterile	Item description	Quant.
<b>Implants – screws</b>			
<b>Cortical screw 4.5 as st</b> 			
2.545.14	S2.545.14	Ø 4.5 mm x 14 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.545.18	S2.545.18	Ø 4.5 mm x 18 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.545.22	S2.545.22	Ø 4.5 mm x 22 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.545.26	S2.545.26	Ø 4.5 mm x 26 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.545.30	S2.545.30	Ø 4.5 mm x 30 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.545.34	S2.545.34	Ø 4.5 mm x 34 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.545.38	S2.545.38	Ø 4.5 mm x 38 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.545.42	S2.545.42	Ø 4.5 mm x 42 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.545.46	S2.545.46	Ø 4.5 mm x 46 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	6
2.545.50	S2.545.50	Ø 4.5 mm x 50 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	3
2.545.54	S2.545.54	Ø 4.5 mm x 54 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	3
<b>Cortical screw 5.5 as st</b> 			
2.549.40	S2.549.40	Ø 5.5 mm x 40 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	2
2.549.45	S2.549.45	Ø 5.5 mm x 45 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	2
2.549.50	S2.549.50	Ø 5.5 mm x 50 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	2
2.549.55	S2.549.55	Ø 5.5 mm x 55 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	4
2.549.60	S2.549.60	Ø 5.5 mm x 60 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	4
2.549.65	S2.549.65	Ø 5.5 mm x 65 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	4
2.549.70	S2.549.70	Ø 5.5 mm x 70 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	4
2.549.75	S2.549.75	Ø 5.5 mm x 75 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	4
2.549.80	S2.549.80	Ø 5.5 mm x 80 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	4
2.549.85	S2.549.85	Ø 5.5 mm x 85 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	4
2.549.90	S2.549.90	Ø 5.5 mm x 90 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5; stainless steel	2

# Set Femur distal groove and notch IH S / IMPLANTS sterile

angle-stable; for screws Ø 4.5/5.5 mm; stainless steel

not sterile packed Set no. 19.600.10 / sterile packed Set no. 19.600.10PS+

Item no.	Item no. sterile	Item description	Quant.
Implants – screws			
Screw head 8.0 for system 4.5			 
8.946.05	S8.946.05	Screw head with conical thread Ø 8.0 mm x 5.0 mm; internal hexagon 3.5; big fragment; stainless steel	5

Instruments			
2.904.09	S2.904.09	Drill bit for quick coupling Ø 3.2 mm; spiral 60 mm; 2-flute; length 195 mm	1
2.904.23	S2.904.23	Drill bit for quick coupling Ø 4.0 mm; spiral 60 mm; 2-flute; length 220 mm	1
2.953.90	–	Gauge with clasp for screws Ø 4.5 mm an Ø 6.0 mm; measuring range 90 mm	1
2.954.01	–	Screw forceps; self holding	1
2.977.02	–	Drill guide for angle-stable screwing; for drills Ø 3.2 mm; length 70 mm; big fragment	2
2.977.15	–	Drill guide for angle-stable screwing; for drills Ø 4.0 mm; length 70 mm	2
2.9406.35	–	Hex screwdriver with handle for screws Ø 4.5 mm – Ø 7.0 mm; wrench width 3.5 mm; length 215 mm	1

Wires			
6.031.16	S6.031.16	Kirschner wire with trocar point and round end; Ø 1.6 mm x 150 mm; stainless steel	5

Container			
19.601.00	–	Perforated autoclavable container with inset for instruments and distal femur plates; Ø 6.0 mm; stainless steel	1

Optional			
2.163.XX	S2.163.XX	Cancellous screw with conical head thread Ø 6.0 mm; fully threaded; self tapping; stainless steel	–
2.941.35	–	Hex screwdriver inset for quick coupling; for screws Ø 4.5 mm – Ø 7.0 mm; wrench width 3.5 mm	–
2.941.35/135	–	Hex screwdriver inset for quick coupling; for screws Ø 4.5 mm – Ø 7.0 mm; length 135 mm; wrench width 3.5 mm	–
2.941.35/195	–	Hex screwdriver inset for quick coupling; for screws Ø 4.5 mm – Ø 7.0 mm; length 195 mm; wrench width 3.5 mm	–
2.9355.35/20	–	Silicone handle with ratchet and AO coupling; length 165 mm; Ø 20 mm	–

#### Applicable with:


- 19.518.00/1 Cable osteosynthesis 1.5 S SET
- 19.604.00 Distal Femur plate 4.5/5.5 as IH aiming device INSTRUMENTS
- 19.045.001 LF tools IH INSTRUMENTS

**Attention:** When using the sterile set (set no. 19.600.10PS+), an additional instrument set (set no. 19.045.001) is required.


# Set Femur distal groove and notch IH S / IMPLANTS sterile

angle-stable; for screws  $\varnothing$  4.5/5.5 mm; stainless steel


not sterile packed Set no. 19.600.10 / sterile packed Set no. 19.600.10PS+



Screws LEVEL 1



Screws LEVEL 2

 Internal hexagon  
Set no. 19.600.10  
not sterile packed