


# Set MF 2.4 tools TOBI® INSTRUMENTS

for screws Ø 2.4 mm; for TOBI®-Drive B8

non-sterile packaging Set no. 19.024.201



## TOBI

Set no. 19.024.201  
non-sterile packaging

Selected items are available individually in non-sterile and sterile packaging.

Item no.	Item no. sterile	Item description	Quant.
<b>Instruments</b>			
2.904.33	S2.904.33	Drill bit with quick coupling Ø 1.8 mm; spiral 30 mm; 2-flute; length 105 mm	1
2.953.24	–	Gauge with clasp for screws Ø 2.4 mm; measuring range 60	1
2.9778.24	–	Drill guide for angle-stable screwing; for drills Ø 1.8 mm; TOBI®-Drive B8; length 46 mm; stainless steel	1
2.9286.08	–	Screwdriver for TOBI®-Drive B8; PPSU handle; length 187 mm; stainless steel	1
10.500.24	–	Combination drill guide with handle for screwing 2.4; for drills Ø 1.8 mm; length 42 mm	1

<b>Wires</b>			
6.031.15	S.6.031.15	Kirschner wire with trocar point and round end; Ø 1.5 mm x 150 mm; stainless steel	1

<b>Container</b>			
19.024.21	–	Container for instruments 2.4 TOBI®-Drive	1

<b>Optional</b>			
10.184.07	–	Guide block for distal radius plate 2.4; 5.08173.XX	–
10.184.09	–	Guide block for distal radius plate 2.4; 5.08193.XX	–
10.185.07	–	Guide block for distal radius plate 2.4; 5.08273.XX	–
10.185.09	–	Guide block for distal radius plate 2.4; 5.08293.XX	–
10.187.07	–	Guide block for distal radius plate 2.4 Watershed; 5.08373.XX	–
10.188.07	–	Guide block for distal radius plate 2.4 Watershed; 5.08473.XX	–

# Set MF 2.4 tools TOBI® INSTRUMENTS


for screws Ø 2.4 mm; for TOBI®-Drive B8

non-sterile packaging Set no. 19.024.201

Item no.	Item no. sterile	Item description	Quant.
Optional			
2.9355.35/20	—	Silicone handle with ratchet and AO coupling; length 165 mm; Ø 20 mm	—
2.9417.B8/75	—	Screwdriver insert; for TOBI®-Drive B8 system; with AO coupling; length 75 mm	—

## Applicable with:

- 19.117.20 Dist. Radius plate 2.4 vas TOBI® T SET
- 19.117.20PS+ Dist. Radius plate 2.4 vas TOBI® T IMPLANTS sterile
- 19.117.30 Dist. radius plate 2.4 Watershed vas TOBI® T Set
- 19.117.30PS+ Dist. Radius plate 2.4 Watershed vas TOBI® T IMPLANTS sterile



**TOBI**

Set no. 19.024.201  
non-sterile packaging