

# Set RoSA® IH T Trochanter stabilization plate / IMPLANTS sterile

Rotation-stable Screw Anchor for extramedullary treatment; variable angle-stable; titanium

non-sterile packaging Set no. 19.612.00 / sterile packaging Set no. 19.612.00PS+



Item no. 7.305.04  
Item no. S7.305.04



Item no. 7.304.04  
Item no. S7.304.04



Item no. 10.269.924



Item no. 10.269.922



All listed implants are available individually in non-sterile and sterile packaging.

Implants in sets, supplied **non-sterile** Set no. 19.612.00

Implants in sets, supplied **sterile** Set no. 19.612.00PS+

 **Internal hexagon**  
Set no. 19.612.00  
non-sterile packaging

| Item no.  | Item no. sterile | Item description  | Quant. |
|---|------------------|---|--------|
| <b>Implants – plates</b>                            |                  |   |        |
| <b>RoSA® trochanter stabilization plate 4.5 vas</b> |                  |   |        |
| 7.304.00  | S7.304.00        | RoSA® trochanter stabilization plate for screws Ø 4.5 mm; 2 head holes; neck length 0 mm; variable angle-stable; titanium | 1      |
| 7.304.04  | S7.304.04        | RoSA® trochanter stabilization plate for screws Ø 4.5 mm; 2 head holes; neck length 4 mm; variable angle-stable; titanium | 1      |
| 7.305.00  | S7.305.00        | RoSA® trochanter stabilization plate for screws Ø 4.5 mm; 4 head holes; neck length 0 mm; variable angle-stable; titanium | 1      |
| 7.305.04  | S7.305.04        | RoSA® trochanter stabilization plate for screws Ø 4.5 mm; 4 head holes; neck length 4 mm; variable angle-stable; titanium | 1      |

| <b>Implantats – screws</b>  |           |   |   |
|---|-----------|---|---|
| <b>Cortical screw 4.5 as st</b>   |           |   |   |
|  |           |   |   |
| 3.545.56  | S3.545.56 | Cortical screw; Ø 4.5 mm x 56 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium | 4 |
| 3.545.60  | S3.545.60 | Cortical screw; Ø 4.5 mm x 60 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium | 4 |
| 3.545.65  | S3.545.65 | Cortical screw; Ø 4.5 mm x 65 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium | 4 |
| 3.545.70  | S3.545.70 | Cortical screw; Ø 4.5 mm x 70 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium | 4 |
| 3.545.75  | S3.545.75 | Cortical screw; Ø 4.5 mm x 75 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium | 4 |
| 3.545.80  | S3.545.80 | Cortical screw; Ø 4.5 mm x 80 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium | 4 |
| 3.545.85  | S3.545.85 | Cortical screw; Ø 4.5 mm x 85 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium | 4 |
| 3.545.90  | S3.545.90 | Cortical screw; Ø 4.5 mm x 90 mm; angle-stable; fully threaded; self tapping; internal hexagon 3.5 mm; titanium | 4 |

# Set RoSA® IH T Trochanter stabilization plate / IMPLANTS sterile

Rotation-stable Screw Anchor for extramedullary treatment; variable angle-stable; titanium

non-sterile packaging Set no. 19.612.00 / sterile packaging Set no. 19.612.00PS+

| Item no.              | Item no. sterile | Item description  | Quant. |
|-----------------------|------------------|---|--------|
| Instruments – Level 2 |                  |   |        |
| 10.269.93             | –                | Fixation screw for RoSA® trochanter stabilization plate; internal hexagon 3.5 mm                      | 2      |
| 10.269.922            | –                | Guide block for RoSA® trochanter stabilization plate; 2 holes   | 1      |
| 10.269.924            | –                | Guide block for RoSA® trochanter stabilization plate; 4 holes   | 1      |
| 10.173.05             | –                | Screw driver for guide block; length 68.5 mm  | 1      |
| 2.9406.25             | –                | Hex screwdriver with handle for screws Ø 2.7 mm – 4.0 mm; wrench width 2.5 mm; length 215 mm; conical | 1      |
| 10.269.94             | –                | Drill guide for variable angle-stable screwing; for drills Ø 3.2 mm for RoSA®                         | 1      |
| 10.269.90             | –                | Bending fixture for RoSA® trochanter stabilization plate  | 1      |
| 10.269.91             | –                | Bending lever for RoSA® trochanter stabilization plate  | 1      |

|           |   |   |   |
|-----------|---|---|---|
| Container |   |   |   |
| 19.611.00 | – | Perforated autoclavable container with inset for instruments / implants RoSA®; trochanter stabilizing plate | 1 |


|           |   |  |   |
|-----------|---|--|---|
| Optional  |   |  |   |
| 2.9406.35 | – | Hex screwdriver with handle for screws Ø 4.5 mm – 7.0 mm; wrench width 3.5 mm; length 215 mm | – |

## Applicable with:

|           |                              |
|-----------|------------------------------|
| 19.606.00 | RoSA® IH T SET               |
| 19.608.00 | RoSA® IH INSTRUMENTS         |
| 19.610.00 | RoSA® IH INSTRUMENTS Removal |

# Set RoSA® IH T Trochanter stabilization plate / IMPLANTS sterile

Rotation-stable Screw Anchor for extramedullary treatment; variable angle-stable; titanium  
non-sterile packaging Set no. 19.612.00 / sterile packaging Set no. 19.612.00PS+



Internal hexagon  
Set no. 19.612.00  
Level 1  
non-sterile packaging



Internal hexagon  
Set no. 19.612.00  
Level 2  
non-sterile packaging